CE
0459

MONT BLANC®

CONTRE-ANGLE D'IMPLANTOLOGIE 20:1
IMPLANTOLOGY CONTRA-ANGLE 20:1
CONTRAANGULO PARA IMPLANTOLOGIA 20:1
CONTRANGOLO D'IMPLANTOLOGIA 20:1
IMPLANTOLOGIE WINKELSTÜCK 20:1
CONTRA-ÂNGULO PARA IMPLANTOLOGIA 20:1

10400X - 10400XL - 10400XLED
10401X - 10401XL - 10401XLED
10403X - 10403XL - 10403XLED

NOTICE D'INSTRUCTION
INSTRUCTIONS FOR USE
NOTA INFORMATIVA DE INSTRUCCIÓN
FOGLIO D'ISTRUZIONI
GEBRAUCHSANWEISUNG
INSTRUÇÕES DE USO

anthogyr

2237, Avenue André Lasquin
74700 Sallanches
FRANCE

Tel +33 (0)4 50 58 02 37
Fax +33 (0)4 50 93 78 60

www.anthogyr.com









Tel SAV: +33 (0)4 50 58 50 53
sav@anthogyr.com

Made in France

NOT 10400- 2013-06

anthogyr

A global solution for dental implantology

SOMMAIRE

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Visual inspection

DM

Medical device
(contra-angle, hand
piece, etc.)

General information

Wear gloves

Fibre optic device for
motor with a light source

REF

Device
reference

Device with light
emitting diode

SN

Device serial
number

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Date of manufacture

Thermo-disinfectable

Name of
manufacturer

Single Use

II.SCOPE OF USE

Surgical treatment for hard organic substances: : drilling, tapping, screwing of dental implants, screwing of the prosthetic parts.

Inappropriate or indirect use could damage this device and be a risk to the user and third parties.

These **MDs** are for professional use only in the dental surgery field.

According to these instructions, the **MD** must only be used by a user with dental medicine experience, for the stated use, and in compliance with current guidelines regarding the prevention of accidents at work and work protection and the instructions in this leaflet. **MDs** must be prepared and maintained only by individuals who have been trained in preventing infections, auto-protection and patient protection.

According to these instructions, users must :

- Only use non-defective work instruments,
- Only use the **MD** in a unit that complies with the guidelines of the EN 60601 standard,
- Follow the correct use instructions,
- Protect oneself and patients and third parties against all dangers,
- Avoid all contamination by the product.

The following situations:

- Inappropriate use,
- Lack of maintenance,
- Use of removable accessories or parts not approved by Anthogyr,
- Using accessories from other devices on this **MD**,
- Change or addition to a **MD** not validated by Anthogyr.

Relieve Anthogyr of all guarantee obligations or other claims.

These **MDs** comply with the Community Directive 93/42/CEE as amended by the 2007/47/CEE directive.

III. GENERAL SAFETY INFORMATION

Prior to use, check that the device has not been damaged in any way and that no parts are missing.

Wear appropriate protection, especially gloves, a mask and glasses.

Do not use the **MD** and inform your distributor or Anthogyr After-Sales (AS) in the following situations :

- visible failure or damage,
- non-existent or insufficient irrigation (at least 50 ml/min, correctly tempered and oriented in an optimum manner),
- scratch or damage to the emitting side of the LED.

Risk of damage to the device and injury, wait until the motor has completely stopped to:

- connect/disconnect the **MD** from the motor,
- turn the assembly/disassembly system on,
- turn the locking/unlocking system for the rotary instrument on,
- handle the rotary instrument.

Assess the risk of septic substances penetrating tissue if used with a central air cooling motor (refer to the motor manufacturer's leaflet).

Warning : this **MD** has a ball bearing system.

The power is maintained :

- For motors having a torque calibration system (like the IMPLANTEO®): you must calibrate the contra-angle as directed by the manufacturer. The motor must have been calibrated for the used contra-angle.
- For any other motors: make sure that the displayed torque values are right.

A tapping or a screwing torque which is too high (above 35N.cm) may increase risks of necrosis or of bone fissuring. Furthermore, the implant and the instrument may be damaged: implant manufacturer's prescriptions have to be respected.

Risk of burning and injury :

- Apply pressure to the **MD** with a rotary instrument (refer to the manufacturer's instructions) without exceeding 2N (about 200g).

Risk of burning related to the LED overheating :
• Avoid all contact between the LED and soft tissue.

Risk of dazzle:

- Do not look at the light output directly (especially on the LED model).
- Do not point the beam directly at the eyes of the patient or third parties.

The lighting level can generally be adjusted in the dental unit.

Risk of electrical discharge due to **MD** being badly connected to an Anthogyr non-compatible system:

In the event of the **MD** being assembled and used in other manufacturer's treatment devices and installations, refer to "Protection against the risk of electrocution", "Leaking current" and "Non-earthing of use part" according to the IEC 60601-1 standard.

Do not use in explosive atmospheres.

IV. TECHNICAL CHARACTERISTICS

IV.1 – Description of MD (Fig.1 and 2)

Code	Description
A	Head
B	Push button
C	Spraying nozzle
D	Collar cartridge
E	Fibre optic or LED outlet
F	Locking system
G	Body
H	Depth stop system
I	Depth stop

IV.2 - Characteristics of MDs

References	10400X	10401X	10403X	10400XL	10401XL	10403XL	10400XLED	10401XLED	10403XLED
Ratio	20 : 1								
Colour code	Green								
Adjustable depth stop system	No	Yes	No	Yes	No	Yes	No	Yes	
Light	No			Yes (light guide)			Yes (LED in the contra- angle)		
Motor connection standard	ISO 3964			ISO 3964 ¹			ISO 3964 ²		
Maximal speed of motor (rpm)	40 000								
Maximum torque of motor (N.cm)	6,5								
Maximum torque of instrument (N.cm)	80 ³								
Rotary instruments according to EN ISO 1797-1 → ø (mm) → Lg max (mm) → ø max in the active part (mm)	Type 1 2,35 40 ⁴ 10 ⁴								
Spray water flow according to ISO 7785-2 (ml/min)	> 50								

S3 10% intermittent mode (1 min on/ 9 min off). The authorised relative load duration is 1 minute per 9 minutes of pause.

IV.3 - Adaptor for LED MD:

Intended to be placed instead of the bulb in dental motors with a INTRAmatic LUX® 1 type connection, the adaptor for LED MD will be supplied with power from the treatment unit.

Electrical characteristics for power of adaptor and/or LED MD:

- Nominal supply voltage for the adaptor for LED MD: 3.4 V CC
- Nominal supply voltage range for the adaptor for LED MD: 3.0 – 3.6 V CC
- Maximum current: 400mA

¹ INTRAmatic LUX® connection system - Kavo® registered trademark

² IMPLANTEO LED® Anthogyr connection system. E-PLUS or E-PLUS Intra specific connection. Can be assembled on a INTRAmatic LUX®¹ type motor for which the bulb has been replaced by an adapter for LED "10510" contra-angles.

³ The torque can be reached providing the motor used is powerful enough: for example, the Implanteo motor by Anthogyr. In general, it is possible to limit the motor torque on the operating unit.

⁴ Indicative values. In the event of the use of longer or bigger rotary tools, the user is responsible for choosing the right operational conditions that avoids all risks to patients or third parties. Follow the instructions for use, in accordance with the tools manufacturer's instructions.

The LED is a semi-conductive element which must only be powered with continuous voltage. In order to ensure adequate operation, it must be properly polarised: the positive terminal of the connection part for powering the LED at the back of the **MD** is closest to the centre (**Fig 3**).

V.INSTALLATION OF MD

The **MD** is supplied non-sterile.

Before first use, the **MD** must be cleaned, decontaminated and sterilised (see § VI)

V.1 – Installation of the adaptor for the LED MD

Check that your treatment unit satisfies the electrical power for the supply voltage (see § IV.3). If the electrical power is not indicated, please contact the manufacturer of your treatment unit or your installer.

- Disconnect the electrical supply for the treatment unit when installing the adaptor for the LED MD. Risk of electric shock.
- Proceed as with replacing the bulb – refer to the motor manufacturer's instructions.
- Insert the adaptor up to the stop.

The polarity of the power for the LED is important: if this does not light up when in use, then use an adaptor for the LED MD with a polarity correction or contact your installer.

V.2 – Connection to the motor

Check that the **MD** is completely dry before connecting it to the motor. Install the MD on the motor connection until ratcheting. For this, keep the motor and the **MD** in the same axle.

or **LED**

: Turn the **MD** until the retractable catch lodges in the E notch of the motor (**Fig 4**).

Apply light traction to the **MD** to check that it is properly attached to the motor connection.

Operational test by switching on the motor; start at low speed, then gradually increase. If you notice overheating, irregularities, vibrations, abnormal noises when using the **MD**, immediately contact your AS.

V.3 – Connection of the irrigation (Fig. 5)

Attach the spraying nozzle to the head of the **MD**.
Connect the extremity of the irrigation pipe to the nozzle entry.

V.4 – Disconnection of motor

Remove the **MD** by keeping it in the motor axle.

In the event of prolonged non use of the device, do not leave it connected to the motor. Risk of damaging the motor.

V.5 – Connection/disconnection of rotary instrument

For all handling of rotary instruments, wear protective gloves. Risk of cutting and contamination.

Only use rotary instruments in good condition. Risk of necrosis.

As the maxi torque given by this **MD** is very high, this may damage the tool's driving part. In any case, the compatibility of your instruments has to be checked. Please follow the manufacturer's instructions and limit the torque of the unit if necessary. We strongly advise that you use instruments having a shaft made out of hardened steel.

Check regularly that the push button is set properly by applying a slight unscrewing torque.

Insertion of rotary instrument

- Apply pressure using your thumb to the centre of the push button and simultaneously insert the rotary tool up to the stop in the claw of the MD.
- Release the thumb pressure on the push button (**Fig 6**)

For tool replacement: check the good condition of the rotary tool with light axial traction.

V.6 – Adjustment of the depth stop

It is imperative that the depth stop is adjusted outside of the patient's mouth.

To adjust the depth stop: move the sliding block backward (1) and move the stop to the required position (2). Loosen the sliding block.

- After each adjustment, check the the stop has been correctly placed using the ruler.
- Before each use, exercise a slight traction on the stop in order to verify that it is correctly installed.
- Do not use the depth stop for tapping: the bone may be damaged in the event of late stop of the driving motor.

The specific depth stops marked AX may only be used for placing Axiom® implants with the Axiom® guided surgery. Verify the position of the depth stop properly in relation to the instrument reference lines.

VI. HYGIENE AND MAINTENANCE

Sterilisation of medical devices must be done by properly trained and protected staff, in compliance with current regulations. The sterilisation protocol must be appropriate to the infectious risk.

Wear appropriate protective clothing: risk of infection and injury.

Only use products for the maintenance of medico-surgical equipment compatible with stainless steel (no chlorinated content).
Forbid antiseptics that are intended for use on skin and mucous.
Forbid products containing aldehyde, alcohol or other products likely to bind proteins.

For each product used: refer to the manufacturer's instructions.
Comply with the concentrations, exposure durations, and life span of products. Do not mix the products and follow the instructions for their disposal.

Spraying nozzles Ref. 10476 (the green ones) cannot be sterilised: change them after each surgery.

VI.1 – Preparation of MD in the place of use

To be done immediately after surgery (maximum time of 2 hrs):

- Disconnect the **MD** from the motor,
- Remove the rotary tool (see § V.6),
- Remove the largest organic contaminants with disinfectant wipes.

VI.2 – Preparation for sterilisation

Disassembly of MD (Fig 1):

- Remove the rotary instrument,
- Disconnect the **MD** from the motor,
- Remove the rotary instrument and the spraying nozzle.
- Pivot the F lever up to the stop (indifferently towards the right or towards the left),
- Pull the A head,
- Remove the collar cartridge D.

VI.2.1 – Manual preparation for sterilisation

Cleaning of irrigation pipe:

Carefully clean the irrigation pipe and remove any impurities and deposits from it using the caustic cleaner spray.

Risk of injury: use protective glasses.

By brushing:

- Brush each part under running water with a soft brush.

LED

Lighting lens (LED or Fibre optic):

- Absolutely avoid all scratches to the LED or the fibre optic,
- Wipe a disinfectant wipe over each of the **MD**'s sub-units,
- Or Spray disinfectant on the **MD** and wipe with a clean cloth.

Or: by Ultrasound

- Use a low frequency ultrasound tank (25 to 50 kHz) and a detergent-disinfectant product compatible with this process and with the **MD**,
- Friction of parts between them or with the tank can cause appearance defects,
- Install light input/output protections on the body (**Fig 7**) to protect the fronts of the fibre optic or the emitting surface of the LED,
- Immerse the components.

**Rinsing and drying:**

- Empty the spray pipe (**Fig 4**).
- Rinse well⁵, then dry each part.

VI.2.2 - Automatic preparation for sterilisation

- Only in a washing and disinfecting machine, appropriate for the treatment of this type of **MD**. Follow the manufacturer's instructions for the device.
- If the device is not equipped with pipes for contra angles / handpieces, then disassemble the **MD** (see § V.2) and immobilise each part. Otherwise, put the **MD** mounted on the pipe.
- The disinfection thermo cycle must be at least 10 minutes at 93°C (203°F).
- Check there is no residue and that all the parts of the **MD** are completely dry at the end of the cycle.
- If the washing and disinfecting machine is near the place of use, it is possible to also do step **VI.1**.

**Re-assembly of MD:**

Only assemble parts with the same serial number⁶.

- Attach the collar cartridge D to A head (only one direction possible) until it completely plates (turn the cog until connection if required),
- Insert the head – collar cartridge unit in the G sleeve (connect a clean tool in the head and turn it slightly so that the cog connects),
- Bring the F lever to the central position, until ratcheting. Make sure that the depth stop system is in the right position.

Systematically check the condition of the front part by applying axial traction.

VI.3 - Lubrication

Follow the spray lubricant manufacturer's instructions.

Wear appropriate protective clothing (splashes, etc.).
Keep away from all sources of heat or ignition.
In particular: do not smoke. Risk of inflammation.

⁵ Temperature < 38°C, distilled water is recommended instead of running water if this has too much chlorine (cf. FD98-135 standard).

⁶ Exception : If the After-Sales Service has to replace a part, the serial number may be different.

After cleaning and before each sterilisation (once a day minimum):

- Remove the rotary instrument,
- Insert the tip to the back of the **MD (Fig 8)**,
- Cover the head of the **MD** with a soft cloth or paper or a wipe.
- Point the head towards the bottom,
- Spray several times by firmly holding the instrument,
- Wipe excess oil with a cloth or a wipe.

VI.4 – Operational test

- Connect the **MD** to a micro-motor, point the head towards the bottom.
- Turn the micro-motor on at low speed for 30s. Gradually bring up to maximum speed.
- Wipe excess oil with a cloth or a wipe

If you notice overheating, irregularities, vibrations, abnormal noises when using the **MD**, immediately contact your AS.

VI.5 - Sterilisation

- The instruments are to be sterilised before first use and after each use.
- Only sterilise cleaned, lubricated and tested instruments.
- We recommend sterilisation that complies with the EN 13060 standard, class B. All other methods are forbidden.
- Independently put each **MD** in a sterilisation sachet that complies with current standards and the size of the **MD**.
- Comply with the space between the sachets and do not overload the autoclave.

• These **MDs** and their accessories (apart from accessories with the logo  which are single use) must be sterilised at 135°C for 18 minutes minimum (sterilisation time).

• Check that the device is completely dry at the end of the cycle.

VI.6 - Storage

Keep the **MD** in a sterilisation sachet in a clean and dry place.

VII.REPAIR

In the event of breakdown, please contact your approved distributor or our after-sales department directly.

Repairs must only be carried out by an approved repairer or by Anthogyr After Sales Department, only with Anthogyr original replacement parts.

For all revisions or repairs, the **MD** must be returned complete and sterile with proof of sterility. It must be accompanied by a document outlining the problem and showing the complete contact details of the user.
The replacement of removable parts is covered for 7 years.

VIII.GUARANTEES

This **MD** is guaranteed parts and labour against all manufacturing defects for **12** months from the date of invoice. This guarantee does not apply to wear.

All changes or additions to the product without the express agreement of Anthogyr render this guarantee null and void.

The guarantee becomes null and void if the technical instructions are not followed.

Anthogyr cannot be held responsible for damage resulting from or which could result from normal wear, use, cleaning or incorrect maintenance, the non-observance of instructions for use or connection, scaling or corrosion, impurities in the water supply system or unusual chemical or electrical influences or non observance of the instructions, maintenance instructions and assembly of Anthogyr and other manufacturer's instructions.

Delivery charges incurred when sending an instrument back to Anthogyr for repair will be paid by the client, even if the repair itself is covered by the guarantee.

Postage and packing fees when returning the instrument to the client are covered by the guarantee.

So that guarantee requests are taken into consideration, please attach a copy of the invoice or a copy of the delivery slip to the **MD**.

IX.ACCESSORIES

To be ordered from your approved distributor.

Description	Reference
Sterilisable Spraying nozzle	DCA 10475
Pack of 10 single use spraying nozzles	DCA 10476D
Sterilisable coupling	DCA 10420
Needle for internal spray	DCA 10430
Clip for irrigation line	DCA 10430
Depth stop holder + depth stop	DCA 10415
Depth stop for Axiom® guided surgery	DCA 10418_SAFE
Ruler	DCA 6458
Protective caps for lighting system	DCA 10436
Anthogyr 500 ml lubricating spray	DCA 1930X
Nozzle for standard ISO type E connection	DCA 1932X
Spray caustic cleaner	DCB 9205
Adaptor for LED/Anthogyr contra-angle	DCA 10510A

X.DISPOSAL

The **MD** must be sterilised before disposal. Risk of contamination of third parties.

Comply with legislation and national standards and guidelines for disposal.

Integrated LED **MD** and adaptor: Comply with legislation, national standards and guidelines for the disposal of small electrical equipment or worn electronic equipment.

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