**MIRAFLUORID®**  
Produktionsgesellschaft mbH  
Schusterstr. 35-23355 Barmstedt/Germany  
T +49 4123 9228-0 F +49 4123 9228-49  
www.dent-a-pharm.de info@dent-a-pharm.de

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**en Instructions for use**

**Patient target group:** Children and adolescents

**Intended users:** Dentist

**Composition:** Qualitative data: Sodium fluoride, bis-(hydroxyethyl)-aminopropyl-N-hydroxyethyloctadecylamine dihydrofluoride, aqueous methacrylate dispersion

Quantitative data: Sodium fluoride, bis-(hydroxyethyl)-aminopropyl-N-hydroxyethyloctadecylamine dihydrofluoride 0,15 % calculated as pure fluoride

**Intended purpose of the product:** Fluoridation of the tooth enamel.

**Indications:** Mirafluorid is a varnish for fluoridation of the tooth enamel and desensitisation of sensitive teeth.

**Contraindications:** Do not use in patients with a known allergy to methacrylates.

**Performance characteristics of the product:** After varnish formation (application of dispersion and dry blowing) on the enamel, the varnish softens and releases fluoride.

**Specific precautions:** None

**Properties:** Mirafluorid is an air-drying dispersion that forms a well-adhering film on enamel and dentin. The adherent film softens when it comes into contact with saliva or water for a longer time. In the process, the fluorides contained in the polymer matrix can diffuse to the dentin or enamel surface and exert their effect there. Due to this behaviour, minimal amounts of fluoride are completely sufficient for optimal fluoridation. Once the varnish is completely soaked, the fluoride release is complete. The varnish matrix then detaches from the enamel, the varnish softens and releases fluoride.

**Application:** Clean surfaces to be treated and, if possible, dry them. Transfer Mirafluorid drop by drop to the brush (it is not necessary to shake up the aqueous dispersion before removal) and touch surfaces gently. Dry the dispersion carefully and extensively, preferably with warm air, until the polymer film has completely formed and is dry. Do not brush teeth for about 24 hours after treatment. When used correctly, the polymer layer will not come off until it is completely soaked and maximum fluoride release has occurred.

**Duration of application:** For caries prophylaxis, treatment should be carried out 2–3 times a year (or 3–4 times if there is an increased susceptibility to caries).

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**Warnings:** Residual risks: The processing instructions and precautionary measures described in the instructions for use of the product must be observed in every instance. If, despite following these instructions, an adverse reaction occurs in the patient or other adverse events occur, please inform us in every case, describing the accompanying circumstances and symptoms as precisely as possible. We follow up on every use on a single patient!

**Residual risks:** The processing instructions and precautionary measures described in the instructions for use of the product must be observed in every instance. If, despite following these instructions, an adverse reaction occurs in the patient or other adverse events occur, please inform us in every case, describing the accompanying circumstances and symptoms as precisely as possible.

**Side effects:** Undesirable side effects of this medical product are extremely rare, provided the product is processed and used properly. If you become aware of any adverse effects – even in doubtful cases – please inform us and describe the accompanying circumstances and symptoms as precisely as possible. We follow up on every indication.

No risks in the event of an overdose are known of.

**Interactions:** No interaction of the product and its metabolites with other medical devices, medicinal products or other substances is known of.

**Storage:** Store Mirafluorid at room temperature and not above 22 °C.

**Minimum shelf life:** Do not use after the expiry date.

**Variants:** [REF 605 810](#)

**Disposal:** Dilute unused liquid with water and dispose of with household waste. Cured residues can also be disposed of with household waste.

**All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient is established.**

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