# anthogyr

A global solution for dental implantology

## SAFE RELAX<sup>®</sup> SAFE REMOVER<sup>®</sup>

DÉPOSE-COURONNE AUTO / MANUEL  
AUTOMATIC / MANUAL CROWN REMOVER  
AUTOMATIK / HANDABZIEHER KRONENABNEHMER  
EXTRACTOR DE CORONAS AUTOMÁTICO / MANUAL  
RIMOVI CORONA AUTOMATICO / MANUALE  
SACA PRÓTESE MANUAL/AUTOMÁTICO

NOTICE D'INSTRUCTION  
INSTRUCTIONS FOR USE  
GEBRAUCHSANWEISUNG  
INSTRUCCIONES DE USO  
FOGLIO D'ISTRUZIONI  
INSTRUÇÕES DE USO

---

## TABLE DES MATIERES

---

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Please wear gloves

### REMARKS

General indications with no danger to people or objects

Please read this!

Visual Examination

### OTHER LOGOS

Autoclave at 135°C

Reference

Serial number

Manufacturer

English

## II. AREA OF APPLICATION

This Medical Device is intended for use in dental treatment in the field of dental surgery, particularly for the fitting of permanent prostheses (single or multiple, bridges). Any other use is prohibited and may be dangerous.

Applied directive:

For the Medical Device , we applied Community Directive 93/42/CEE, amended 2007/47 transposed by the French regulations on medical equipment (Official Bulletin - L169 - 12th of July, 1993).

In accordance with these provisions, the Medical Device must only be used by someone with experience of dental medicine for the application described, and in compliance with the provisions in effect in relation to the prevention of accidents at work, protection of employment and the instructions for use. This medical device must be prepared and maintained only by people trained in the prevention of infection, self-protection and protection of patients.

In accordance with these provisions, the user must :

- only use working instruments without defects,
- only use the the device on drive units which comply with the directives of Standard : EN 60 601 / 93/42/CEE,
- observe correct use,
- protect themselves and the patient or third party against any danger,
- prevent all contamination by the product.

### III. GENERAL SAFETY INFORMATION

Before use, check that the Medical Device is not damaged and that no part of it is missing.

- Use protective gloves, glasses and a mask.
- When inserting or changing the instrument, make sure that the Medical Device is stationary.
- After insertion of a instrument, check that it is correctly held by a slight rotary movement.
- Do not handle the instrument during operation.
- Only operate the locking catch when the motor is stopped.
- The attachments should be connected or disconnected and away from the patient's mouth.
- The Medical Device should be connected or disconnected to the motor when the motor is stationary.

In the event of visible malfunction or damage, immediately stop using the instrument and inform your approved distributor or the manufacturer.

In the event of any questions on the equipment, contact either your approved distributor or the manufacturer.

No alteration or addition is to be made to the product without Anthogyr's express agreement.

Only use accessories designed for this device. Do not use the accessories of Medical Device on another type or another brand of Medical Device.

### IV. TECHNICAL SPECIFICATIONS

#### 4.1 AUTOMATIC CROWN REMOVER :

| REFERENCE   | DCB 6950C          | DCB 6961C |
|---|--------------------|-----------|
| Weight (g)  | 260                | 288       |
| Length MINI / MAXI (mm) of Crown-remover (excluding tool) | 132 / 140          |           |
| Maximum diameter (mm)                                     | Ø 20               |           |
| Motor connection standard                                 | EN 23964           |           |
| Minimum recommended motor speed (rpm)                     | 5 000              |           |
| Maximum recommended motor speed (rpm)                     | 25 000             |           |
| Minimum stroke frequency at 5000 rpm                      | 200 Strokes/min    |           |
| Maximum stroke frequency at 25000 rpm                     | 1000 Strokes/min   |           |
| Number of intensity adjustment settings                   | 5                  |           |
| Number of hooks   | 4                  |           |
| Number of wire tips                                       | /                  | 3         |
| Attachments to use  | Only ANTHOGRYR !   |           |
| Attachment connection                                     | Licensee ANTHOGRYR |           |

#### 4.2 MANUAL CROWN REMOVER :

| REFERENCE  | DCB 1820A          | DCB ATD310B |
|--|--------------------|-------------|
| Weight (g)   | 177                | 279         |
| Maximum length MAXI (mm) of Crown-remover (excluding tool) | 190                |             |
| Number of hooks  | 4                  |             |
| Number of wire tips  | /                  | 3           |
| Attachments to use   | Only ANTHOGRYR !   |             |
| Attachment connection                                      | Licensee ANTHOGRYR |             |


## V. PRODUCT DESCRIPTION

### 5.1 AUTOMATIC CROWN REMOVER :

- 1 Motor connection
- 2 Intensity adjustment band
- 3 Attachment connection (Bayonet)
- 4 Handle

### 5.2 MANUAL CROWN REMOVER :

- 1 Attachment connection (Bayonet)
- 2 Rod
- 3 Counterweight

English

## VI. USING SAFE RELAX

The Medical Device is supplied **unlubricated and unsterilised**. Before first use, the Medical Device should be **cleaned, disinfected, lubricated and sterilised**. (Cf §7: "Hygiene and maintenance")


Before use, check that The Medical Device is not damaged and no part is missing.

### 6.1 Connecting to the motor :

(Automatic Crown-remover only)

**WARNING:** Medical Device should be connected to the motor only when the motor is stationary and away from the patient's mouth.

Lock the Medical Device to the motor coupling (as for EN 23964). To do so, keep the motor and The Medical Device on the same axis.

For certain motors, it is necessary to operate the side button to allow connection.


Pull The Medical Device slightly before use in order to check that it is securely connected to the motor coupling.

## • Operating test

Activate the motor in the normal direction (the device will not function if the motor is in "reverse" mode), to engage the system, exert slight traction on an attachment which has already been affixed (see § 6.3 - Connecting an attachment) until you feel micro-strokes.

Check to make sure that the Device is producing regular strokes.


NOTE: If you notice heating, irregularities, vibrations or abnormal sounds during the operation of Medical Device, immediately contact your approved technical department.

## 6.2 Disconnecting from the motor :

WARNING: Medical Device, should only be disconnected when the motor has stopped and away from the patient's mouth.

Remove the Medical Device whilst keeping it in the motor axis.

For certain motors, it is necessary to operate the side button to disconnect the Medical Device.

NOTE: In the event of prolonged non-use of Medical Device, do not leave it connected to the motor so as to avoid oil leakage into the motor. The motor could be damaged.

## 6.3 Adjustment :

(Automatic Crown-remover only)

### Stroke frequency adjustment:

(Automatic Crown-remover only)

The stroke frequency is directly linked to the rotation speed of the motor. Maximum speed: 25000 rpm. The higher the input speed, the higher the stroke frequency will be (close strokes). Conversely, a decrease in motor speed will produce widely spaced strokes.

The table below illustrates the relationship between motor speed and stroke frequency:


| rpm         | 10000 | 12000 | 15000 | 20000 | 25000 |
|-------------|-------|-------|-------|-------|-------|
| Strokes/min | 400   | 500   | 600   | 800   | 1000  |

### Stroke intensity adjustment:

(Automatic Crown-remover only)

The Automatic Crown-remover is fitted with a notched adjustment band which has 5 intensity settings, making it adaptable to the requirements of the specific clinical case.

The settings are identified as follows, from the lowest to the highest:

NB: The notch for the desired setting (5 shown here) should be aligned with the "Safe Relax" logo.

After the clinical case has been analysed, it is recommended that a low frequency and intensity be used to initiate treatment, and gradually increase both settings if necessary. Before any operation is carried out, it is strongly recommended to verify that the medical device is working properly, away from the patient's mouth.

#### **6.4 Connection/disconnection of the instrument :**

(Applies for both Automatic and Manual Crown-remover)

It is preferable to wear protective gloves for all tool handling. Cut risk.

Check the condition of the rotary instruments used and handle them cautiously and carefully.

It is essential that the motor is stopped when inserting or extracting the instrument and away from the patient's mouth (except wires).

#### **Inserting and locking the instrument:**

- Insert the attachment (Hook or ATD Wire-Tip) into the bayonet by



pushing it fully into the recess.

- Give a quarter-turn to the left to lock the spigot into its ratchet position.
- Check that the instrument is properly attached, through a slight axial movement with each tool change.

## 6.5 Use of the hooks:

(Applies for both Automatic and Manual Crown-remover)

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3 Flat hook with non skidding :  
Réf 5929

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## 6.6 Use of ATD Wire-Tips:

(Applies for both Automatic and Manual Crown-remover)



### Bridges :

1 • Insert the metal wire in the interdental space, as near as possible to the crown, either from palatal or labial side.


Close the loop by inseting the sharp end in its groove.

- 50 mm wire for the anterior.
- 60 mm wire for the posterior.


2 • Attach the Wire-Tip to the Crown-Remover.  
(See §6.4 Connecting the attachment)

3 • Activate the Crown-remover.


4 • To remove a long span bridge :  
Fasten several wire-tips to avoid moving them.  
Move the crown-remover from one wire to the next in the given order to divide the traction power.


#### Inlays-cores :

• To remove an artificial tooth stump (inlay-core),  
bore an horizontal labiolingual hole of a least 1 mm  
diameter.  
Insert the metal wire and proceed as with a  
bridge.

#### Richmond :

• To remove a Richmond proceed in the same  
way.

#### CERAM wire tip :


To remove a small ceramic bridge (without  
metal frame) or a ceramo metallic bridge  
or a temporarily fixed bridge on parallel  
abutments :

Take the ceramic-tip with a 95 mm screwed  
wire for a 3 el. bridge, or with a 150 mm  
screwed wire for a 4 el. bridge.


1 • Insert the metal wire in vestibulo-palatin  
or lingual way in the first interdental space,  
starting with the most anterior tooth of the  
bridge.


2 • Place the wire in anterior/molar direction  
through the hole of the tip.


3 • Insert the wire in the following interdental  
space in the same vestibulo-palatin or lingual  
direction.

4 • Close the metal loop by passing the sharp  
end in its groove.

5 • Do it a third time if necessary.  
Fasten the adaptor on the wire-tip and activate  
the crown-remover.

#### To change the wire :

• Unscrew the screw preferably (above a plate to avoid losing  
it).

• Place the new wire short and unsharp side in the threading  
spot (marqued 0). Screw thoroughly.

## VII. HYGIENE AND MAINTENANCE

Re-sterilization of reused medical appliances has to be done by someone who has been trained and is protected, and the regulations in force have to be adhered to. The re-sterilization protocol must be adapted to the infectious risk.

### **GENERAL RULES**

- Refer to the manufacturer's instructions for each product used. Take particular care to adhere to the instructions regarding concentrations, lengths of exposure, the replacement of solutions, and the lifespan of the products.
- Only use products intended for the maintenance of medical and surgical equipment.
- Do not use antiseptics which are designed to be used on the skin and mucous membranes. Use only products which are pH neutral or alkaline (between 2.5 and 9.0 inclusive). Do not use products containing Aldehydes, alcohol or other products liable to set proteins.
- Use only products compatible with stainless steel (chlorine-free ingredients).
- Wear suitable protective clothing. In order to avoid any risk of infection and injury, it is vital to wear protective gloves. Furthermore, the motor must always be stopped before handling the equipment.
- Never immerse the Medical Device in any kind of solution. Never clean the Medical Device in an ultrasonic tank. The motor must always be disconnected or stopped before handling the equipment in any way.
- Sterilization of reused medical appliances has to be done by someone who has been trained and is protected, and the regulations in force have to be adhered to. The sterilisation protocol must be adapted to the infectious risk.
- Never combine products. Adhere to the guidelines regarding the disposal of used products.
- The following procedures are to be carried out after each operation: Pre-disinfection, cleaning, disinfection, lubrication and sterilisation.

### **7.1 Pre-disinfection :**

Pre-disinfection has to be made immediately after the operation.

- External disinfection using a spray or microbiologically-controlled disinfectant wipes.

Never immerse the Medical Device in any kind of solution. Never clean the Medical Device in an ultrasonic tank.

### **7.2 Cleaning :**

- Rinse the device in demineralised water.
- Carefully dry each part.

### **7.3 Disinfection :**

- With disinfectant wipes (recommended) or a spray. Follow the manufacturer's instructions for use.

### **7.4 Lubrification :**

- This should be done before first use.
- It also must be done before every sterilization.

- Keep away from any heat or fire source. Do not smoke!
- Wear a protective mask.
- Remove the instrument
- Insert the spray into the Dispositif Médical.
- Make sure that the nozzle is the right one.
- Cover the head of the Medical Device with a small decontaminating cloth.
- Put the head upside down.
- Spray several times until the liquid coming out of the Medical Device is clear.
- Wipe the oil excess off with a rag or with a small cloth.

## 7.5 Operating test :

This test must be performed prior to each sterilisation.

Connect the Medical Device to the motor ; point the head of the Medical Device downwards. Operate the Medical Device slowly (5 000 rpm motor) for 30 seconds in order to get rid of the excess oil, then at full speed (25 000 rpm motor) for 30 seconds. Wipe the Medical Device with a disinfectant wipe if the oil has dripped.

NOTE : if you notice a rise in temperature, irregularities, vibrations, or abnormal noises during operation of the Medical Device, contact your approved technical service immediately.

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The Medical Device is supplied unlubricated and unsterilised. The Medical Device should be sterilised before they are first used and following each use.

- Only sterilise instruments which have previously been disinfected, cleaned, lubricated and tested.
- The Device should be sterilised at 135°C at 2.13 bars (275°F at 30.88 psi), in a steam steriliser only, for a minimum duration of 20 minutes (time for which sterilisation must be continued).
- We strongly recommend the use class B autoclave .

- Any other method of sterilisation should be avoided.
- Read the instruction leaflet provided by the steriliser manufacturer.

- Remove the instrument from the Medical Device prior to sterilisation.
- Use sterilisation pouches suitable for Medical Device and the steriliser, in accordance with the standard NF EN 868. Always use one single Medical Device per pouch.
- Adhere to the space specified between pouches and do not overload the steriliser.
- Make sure that the device does not have any areas of corrosion or cracks, and check that it is operating properly. Ensure that the product is dry ; if necessary, dry any residual water with medical quality pressurised air.
- After each sterilisation cycle, check that there is no water remaining on the inside and outside of the packaging. Make sure that the flow indicator has changed to the correct color.
- In order to avoid any retention of water, place the pouch in the steriliser in such a way that any concave parts are face down.
- Keep the devices in the sterilisation pouches away from light, moisture and any contamination. Follow the manufacturer's recommendations as seen on the packaging.
- The duration for which the device is kept after sterilisation should not exceed 1 month. Label the devices, specifying the expiration date.

After the expiration date, repeat the cleaning and sterilisation procedure.

## VIII. REPAIR

In the event of breakdown, please contact your approved distributor or the Anthogyr after sales service department directly.

All repairs must be carried out with parts and subassemblies approved by the manufacturer.

Repairs must only be carried out by an approved distributor or by the after sales service department at the factory.

For any revision or repair the instrument must be returned complete and sterile with proof of sterility.

It must be sent together with a document stating the problem as well as the complete name and address of the dentist.

For any claims to be considered under the warranty, please enclose a copy of the invoice or the delivery note with the equipment.

Replacement of parts is guaranteed for 7 years after the product ceases to be marketed.

## **AFTER-SALES DEPARTEMENT CONTACT DETAILS**

Service S.A.V.

ANTHOGYR

2237, Avenue André Lasquin – 74700 Sallanches - FRANCE

Direct line : +33 (0)4 50 58 50 53

Mail : s.serra.sav@anthogyr.com

## **IX.ACCESSORIES**

To be ordered from your approved distributor.  
Use ANTHOGYR attachments ONLY.

| <b><u>Description</u></b>         | <b><u>Reference</u></b> |
|-----------------------------------|-------------------------|
| Plain ATD wire-tip                | ATD012                  |
| CERAM ATD tip                     | ATD022                  |
| ATD cables, 50 and 60 mm (par 4)  | ATD014                  |
| ATD cables, 95 and 150 mm (par 4) | ATD024                  |
| M2x3 screw for wire-tip           | ATD080                  |
| Screwdriver for wire-tip          | ATD082                  |
| Set of 4 cables 50 mm long        | ATD014-50               |
| Set of 4 cables 60 mm long        | ATD014-60               |
| Bridge hook                       | 5925                    |
| Plain hook for anterior teeth     | 5927                    |
| Elbow hook for posterior teeth    | 5928                    |
| Flat hook with non-slip grooves   | 5929                    |
| Disassembly key for device head   | 87M                     |
| ANTHOGYR lubricating spray        | 1930X                   |
| ISO type E connection tip         | 1932X                   |

## **X. WARRANTY**

This **MD** is guaranteed parts and labour against all manufacturing defects for **12** months from the date of invoice.

This guarantee does not apply to wear and tear parts.

All changes or additions to the product without the express agreement of ANTHOGYR render this guarantee null and void.

The guarantee becomes null and void if the technical instructions are not followed.

ANTHOGYR cannot be held responsible for damage resulting from or which could result from normal wear, use, cleaning or incorrect maintenance, the non-observance of instructions for use or connection, scaling or corrosion, impurities in the water supply system or unusual chemical or electrical influences or non observance of the instructions, maintenance instructions and assembly of ANTHOGYR and other manufacturer's instructions.

Delivery charges incurred when sending an instrument back to ANTHOGYR for repair will be paid by the client, even if the repair itself is covered by the guarantee.


Postage and packing fees when returning the instrument to the client are covered by the guarantee.

So that guarantee requests are taken into consideration, please attach a copy of the invoice or a copy of the delivery slip to the **MD**.

## **XI. CONDITIONS OF STORAGE AND TRANSPORT**

---

## **XII. DISPOSAL OF THE PRODUCT**

---

As far as is currently known, the product does not contain any substances which are harmful to the environment.

The product must be sterilised before being disposed of.

Observe national rules with regard to disposal.

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