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## Instructions for use

36% phosphoric acid wellsamed etching gel with optimum thixotropic flow characteristics. It serves as an etching of tooth substance for enamel and total etch technique. Taking the safety and application instructions into consideration, there are no known side effects to date. Could remove underfilling material by prolonged usage of etching gel.

### Application:

Filling up of empty applicator syringe (only applicable for wellsamed Etching Gel Maxi):

The empty applicator syringe should be locked onto the adapter of the supply syringe (180° clockwise). Fill it up by carefully and slowly pressing the plunger. Eventually pull simultaneously the plunger and the syringe, paying attention not to press the plunger off the syringe. Lastly, turn the adapter to remove it from the applicator syringe and the 50 ml syringe should be tightly closed. Fix the canula onto the applicator syringe.

### Preparation of the applicator syringe:

Take the applicator syringe out of the packaging and remove the cap. An application canula should be screwed onto the syringe (clockwise). It is now ready to be used by carefully pressing the syringe.

### Preparation:

The area of the tooth substance to be etched should be clearly isolated.

- 1) The use of a coffer dam is recommended. Any contact with soft tissues (skin, oral mucosa, eyes, etc.) is to be avoided. In case of contact with soft tissues, immediately rinse with plenty of water and seek medical advice.
- 2) Before treatment, clean the tooth with fluoride- and grease-free polishing paste with the help of a rubber cup.
- 3) Remove cleaning residues with water and dry the surface with air pressure.
- 4) Prepare the respective cavity, chamfer edges to achieve permanent filling edges. Rinse the cavity with water and dry with air pressure. Avoid overdrying.
- 5) Protect the dental pulp in the cavity (e.g. cover with a thin layer of underfilling). The cavity should be dry and free from impurities but avoid overdrying the dentin.

### Doctor's process:

- 6) Apply the wellsamed etching gel on enamel edges of the cavity and when necessary on the dentin:

**Recommended resting time:**  
15-30 sec. on the enamel  
10-15 sec. on the dentin  
50-80 sec. on milk tooth and tooth with high fluoride content

- 7) Minimum 20 sec. rinsing with water.
- 8) Dry with oil- and water-free air pressure. The etched surface is not to be disturbed until the next process or contaminated by saliva.  
In case of contamination, repeat the etching process (approx. 10 sec.)

**The etched enamel should be chalk white!**

After etching it can be proceeded applying adhesives etc., according to the manufacturer's instructions for use.

### SAFETY INSTRUCTIONS:

- R 34 Causes burns.  
S 25 Avoid contact with eyes.  
S 26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.  
S 37/39 Wear suitable gloves and eye/face protection.  
S 60 This material and its container must be disposed of as hazardous waste.

**Content: 36% ortho-phosphoric acid.**

- Only for use in dentistry!
- Not to be used after expiry date.
- Syringe must be immediately closed after use.
- Keep out of reach of children.
- Store at max. temperature of 25°C.

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