**ViscoStat[®], ViscoStat[®]
Wintermint,
Astringent[®] &
Astringent[®] X
Coagulative
Hemostatic Agents****Descriptions**

ViscoStat and ViscoStat Wintermint hemostatics are 20% ferric sulfate equivalent solutions with inert binding agents in a viscous, aqueous carrier. It is kind to hard and soft tissue with a pH of approximately 1.0. ViscoStat Wintermint has a mild mint flavor.

Astringent hemostatic is a 15.5% ferric sulfate equivalent solution in an aqueous carrier with a pH of approximately 1.0.

Astringent X hemostatic is a 12.7% iron ion aqueous solution containing ferric subsulfate and ferric sulfate with a pH of approximately 1.0.

Indications for Use

For use in arresting various degrees of oral capillary bleeding or sulcular fluids during dental treatment and oral surgery including: fixed prosthodontics, restorative-operative, and periodontal treatment. Astringent or ViscoStat may also be used for retrofills, tooth impactions, gingivectomies and as a "fixative" for pulpotomies.

Astringent X is used in situations requiring more profound hemostasis and is not as gentle to the tissues as the ViscoStat and regular Astringent yet is gentler than other iron sulfate solutions. NOTE: Diluted Astringent X does not equal ViscoStat or Astringent.

Directions for Use - Impressions

1. For 1.2ml syringe
 - a. Remove Luer cap.
 - b. Securely attached working tip of choice (Metal Dento-Infusor or Blue Mini Dento-Infusor tip).
 - c. Verify flow prior to applying intraorally.
2. For IndiSpense Syringe
 - a. Remove Luer cap from Indispense syringe.
 - b. Attach a 1.2ml syringe to the male threads of the IndiSpense syringe.
 - c. Depress IndiSpense plunger while guiding 1.2ml syringe plunger to desired fill.
 - d. Separate syringes and re-cap IndiSpense syringe.
 - e. Securely attach working tip (Metal Dento-Infusor or Blue Mini Dento-Infusor tip).
 - f. Verify flow prior to applying intraorally.

NOTE: To avoid cross contamination, use new syringes and tips for additional volumes.


3. Using a palm grasp, slowly express solution while rubbing firmly against the cut tissue. (Fig. 1)
4. Continue around the full circumference of the preparation, expressing and rubbing solution into cut tissue surface to stop bleeding and remove coagulum build up.
5. As excess dark coagulum appears, an assistant should express a gentle stream of water over the area and suction. Hemostasis is achieved when new coagulum stops forming.
6. When hemostasis is obtained, use a firm air/water spray to clean preparation and check for hemostasis. If bleeding occurs, repeat step 3-5 above to bleeding area and re-check with air/water spray.

NOTE: Recommended contact time for Ultradent ferric sulfate hemostatics is 1-3 minutes and up to a maximum of 10-20 minutes.

7. Displace tissue by packing size appropriate Ultrapak or displacement cord into sulcus. For optimum displacement, place a small amount of hemostatic solution (Astringent, ViscoStat, and aluminum chloride or aluminum chloride solutions) in dappen dish and soak cord prior to packing.

NOTE: Do not use "Ultrapak E" or other epinephrine preparations with any ferric sulfate hemostatic including ViscoStat, ViscoStat Wintermint, Astringent, Astringent X as a blue/black precipitate will occur.

8. Thoroughly rinse preparation and surrounding tissue. Wait 1-3 minutes before removing cord.

(Fig. 1)

NOTE: Ideally, the preparation should be cleaned with a scouring product such as Consepsis Scrub before making the final impression.

9. Remove cord, rinse with a firm air/water spray, check for hemostasis, and make impression.

Directions for Use – Direct Bonded Restorations

1. Follow Steps 1 and 2 in the Impressions section to prepare the syringe.

NOTE: If adequate space at bleeding site allows, follow Steps 3-9 in "Impression" section above.

2. Soak Ultrapak or cord in hemostatic solution.

NOTE: Do not use "Ultrapak E" or other epinephrine preparations with any ferric sulfate hemostatic including ViscoStat, ViscoStat Wintermint, Astringedent, Astringedent X as a blue/black precipitate will occur.

3. Displace tissue and control sulcular fluid by gently packing size appropriate cord into sulcus.

4. Thoroughly rinse preparation and surrounding tissue using a firm air/water spray to clean and check for hemostasis (Read CLEANSING NOTE below).

5. Scrub preparation surface with pumice or Ultradent's Consepsis Scrub.

6. Wait 1-3 minutes before removing cord.

NOTE: Cord may be removed or left in place to protect soft tissue.

7. Rinse again with a firm air/water spray and dry. Recommended contact time for Ultradent ferric sulfate hemostatics is 1-3 minutes and up to a maximum of 10-20 minutes.

CLEANSING NOTE: When the tooth is not thoroughly cleaned, residual hemostatic agent or coagulum left on the tooth surface or surrounding tissue may contaminate dentin and/or enamel substrate jeopardizing the bond and seal causing microleakage. Hemosiderin from residual blood in the sulcus or in the coagulum can move between the restoration and tooth producing a dark stain on the underlying preparation. This may occur within days or weeks of placing the restoration and will require retreatment and replacements. When using a total etch phosphoric acid bonding system, the phosphoric acid with help in cleaning the residual hemostatic from the surface however cleaning with a pumice or Consepsis Scrub is always recommended. When using a self-etch bonding system, the preparation should always be scoured with pumice or Ultradent's Consepsis Scrub using a rubber cup or STARbrush.

8. Apply bonding agent and restorative as per manufacturer's instructions.

Directions for Use – Indirect Bonded/Luted Restorations

1. Follow Steps 1-7 for the "Direct Bonded Restorations." Read and understand "CLEANSING NOTE."

2. Once preparation is thoroughly clean, proceed with cementing the final restoration.

NOTE: Temporary cements and/or protein mucins from saliva can contaminate preparations. When in mineral-based hemostatics, (i.e. - Ultradent's ferric sulfate hemostatics) come in contact with protein mucins, they are more firmly attached to the preparation requiring a more aggressive scour-type cleansing after hemostasis is obtained. To avoid poor bond strengths and microleakage as explained above in "CLEANSING NOTE", cleaning with Ultradent's Consepsis Scrub and a rubber cup or STARbrush should be done prior to taking the final impression and prior to final cementation.

If dark stain is present on preparations upon removal of the provisional restoration, it is most commonly due to non-sealing provisional cement allowing the hemosiderin to migrate between the preparation and the provisional restoration. This may occur within days or weeks of placing the provisional restoration. This stain can be very tenacious and requires firm scraping and/or ultrasonic instrumentation followed by an aggressive scouring with pumice or Consepsis Scrub, to remove. This problem can be avoided by using a hydrophilic provisional sealing cement such as UltraTemp regular.

Directions for Use – Vital Pulpotomies (Astringedent or ViscoStat ONLY)

1. Follow Steps 1 and 2 in the Impressions section to prepare the syringe.

2. Scrub firmly at site of bleeding. As excess dark coagulum appears, an assistant should express a gentle stream of water over the area and suction.

3. Remove all hemostatic and extraneous coagulum using an excavator, and firm air/water wash. Recommended contact time for Ultradent ferric sulfate hemostatics is 1-3 minutes and up to a maximum of 10-20 minutes.

4. Apply a thin layer of ZOE, followed by a thin layer of Ultra-Blend Plus keeping both to a minimal thickness so dentin surface is available for bonding.

5. Apply bonding agent and restorative as per manufacturer's instructions.

Directions for Use – Endodontic Surgical Site (Astringedent or Viscostat ONLY)

1. Follow Steps 1 and 2 in the “Impressions” section to prepare the syringe.
2. Using a palm grasp, slowly express solution while rubbing firmly against the bleeding bone and adjacent bleeding areas. As excess dark coagulum appears, an assistant should continuously rinse with sterile water over the area and suction. Recommended contact time for Ultradent ferric sulfate hemostatics is 1-3 minutes and up to a maximum of 10-20 minutes.
3. Once hemostasis is obtained, place retrofilling.
4. If bone grafting product is used, carefully remove all hemostatic and extraneous coagulum from the boney crypt using a curette and sterile water to re-establish bleeding for the healing process. Small quantities of coagulum within the capillaries and small vessels are not problematic.
5. Close surgical site.

Precautions and warnings:

1. Carefully read and understand all instructions, precautions and warnings before use.
2. Recommended contact time for Ultradent ferric sulfate hemostatics is 1-3 minutes and up to a maximum of 10-20 minutes.
3. Hemostatics and/or blood hemostatic mixtures can prevent quality adhesion and lead to staining microleakage under “direct bonded” restorations. They must be thoroughly cleaned using a firm air/water spray and pumice or Consepsis scrub. Residual mucins, coagulum, or residual hemostatics will contaminate the bonding interface and/or prevent polymerization/set of resins.
4. Temporary cements and protein mucins can contaminate preparation surfaces for indirect “bonded/luted” restorations. Mineral-based astringents such as ferric sulfate, can cause these materials to become more firmly attached to the tooth surface. Scouring the surface with Consepsis Scrub or pumice is imperative to ensure the preparation site is absolutely clean. This will ensure the bonding is not compromised which reduces the possibility of microleakage and staining under the restoration.
5. For definitive bonding and luting of restorations, thoroughly scouring and cleaning the preparation is crucial to avoid microleakage and decrease in bond strength.
6. Hemostatics must be thoroughly washed from the preparation site with a firm air/water spray to avoid reaction with polyester materials and thereby compromising the surface set of the impression.
7. When using self-etch bonding agents, the tooth/preparation surface must be scoured with pumice or Consepsis Scrub and thoroughly washed before application. This is not necessary when using a phosphoric etch bonding system or when using conventional glass ionomer, zinc phosphate, or similar cements.
8. To avoid blue discoloration, these ferric sulfate hemostatics should not be mixed with epinephrine or any other epinephrine-impregnated displacement cord.
9. Always thoroughly wash hemostatics from tissues before applying other chemistries such as peroxide-based products.
10. Do not use on patients allergic to ferric sulfate.
11. Ultradent’s hemostatic agents are designed for intraoral use. Verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check. Use only recommended tips.
12. Soft tissue discoloration is temporary.
13. All syringe tips and empty syringes are disposable product and for single use only to avoid cross-contamination. Fill an empty syringe with the amount of material needed for the individual patient. Dispose of syringe after use.
14. Prefilled syringes can be used several times, when protected during each use by Syringe Covers. Please note the instructions for use of the Syringe Covers. Re-cap syringe with the Luer Lock Cap and disinfect syringe with an intermediate level.
15. Do not allow product to be ingested.
16. Avoid contact with clothing.
17. Use protective clothing and eye shield when loading and handling these ferric sulfate hemostatics.
18. Keep out of reach of children.
19. Dispose of used tips and empty syringes properly.

NOTE: For MSDS and additional information about using ViscoStat, ViscoStat Wintermint, Astringedent, Astringedent X, or related products, please go to www.ultradent.com

**ViscoStat® , ViscoStat®
Wintermint,
Astringedent® &
Astringedent® X****Koagulierende, blutstillende
Mittel****Beschreibungen**

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



Ultradent syringes have an expiration date stamped on the side of the syringe consisting of one letter and three numbers. The letter is a lot number used for manufacturing purposes and the three numbers are the expiration date. The first two numbers are the month, and the third number is the last number of the year.

Viscostat		Astringedent		Astringedent X		Ultrapak Cord						
Flammability		Flammability		Flammability		Flammability						
Health		Reactivity	Health		Reactivity	Health		Reactivity	Health		Reactivity	
For product MSDS please see our website: www.ultradent.com				For product MSDS please see our website: www.ultradent.com								
HAZARD RATING												
4 = Severe												
3 = Serious												
2 = Moderate												
1 = Slight												
0 = Minimal												

Manufactured by:
Ultradent Products Inc.
505 West 10200 South
South Jordan, Utah 84095 USA

Made in USA

For professional use only. Keep out of reach of children.

For immediate reorder and/or complete descriptions of Ultradent's product line, refer to Ultradent's catalog or call Toll Free **1-800-552-5512**. Outside U.S. call (801) 572-4200, or visit www.ultradent.com.

For Canadian Import Only
Clinical Research Dental Supplies
London ON N6A 1M6
Viscostat / Viscostat Wintermint - NPN 80016405
Astringedent X - NPN 80006172
Astringedent - NPN 80024496

 **EC REP** UP Dental GmbH
Am Westhover Berg 30
51149 Cologne Germany

#71276.13 022513

 **ULTRADENT**
PRODUCTS, INC.

© Copyright 2011 Ultradent Products, Inc. All Rights Reserved.