CHARISMA®

CHARISMA®
6603406405

- DE Gebrauchsanweisung
- GB Instructions for use
- FR Mode d'emploi
- ES Instrucciones de uso
- IT Istruzioni per l'uso
- PT Instruções de uso
- NL Gebruksaanwijzing
- SE Bruksanvisning
- DK Brugervejledning

- NO Bruksanvisning
- FI Käyttöohjeet
- GR Οδηγίες Χρήσης
- CZ Návod k použití
- HU Használati utasítás
- LV Lietošanas instrukcija
- LT Naudojimo instrukcija
- PL Instrukcja obsługi

Heraeus

CHARISMA®

Universelles lichthärtendes Microglass® Komposit, röntgenopak

Farben:

A1, A2, A3, A3.5, A4, A1, B1, B2, B3, C2, C3, C4, D3, O2, O4, O43, O43.5, O82, Incisal, SL, SLO

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CHARISMA® is a lightning, radio-opaque Microglass® composite. Due to the grain size of the glass filler the final polishing may easily be realized and in connection with adequate bonding systems it is meant for adhesive filling therapy.

Composition:

CHARISMA® is based on a BIS-GMA matrix and contains 58% filler (by volume), which is:

• Bariumaluminum fluoride glass (0,02–2 microns)

• Highly disperse Siliciumdioxide (0,02–0,07 microns)

Application:

• Filings of the classes I–V according to Black

• Adjustments of shapes and shades

• Splitting of teeth which are traumatically loose

• Diastema closure

• Restorations of primary teeth

• Adjustments of developmental teeth defects (e.g. enamel dysplasia)

Kind of application:

In order to prepare the treatment it is recommended to clean the hard dental substance by means of a polishing paste which is free of fluoride. Prior to starting the treatment the shades should be selected by means of the attached colour scale as long as the teeth are still humid. The Heraeus Kulzer Plasmacat instrument as well as devices consisting of Wironit or plastic are specifically suitable for the application.

1. Preparation of the cavity according to the general regulations of composite filling therapy. Then clean and dry the cavity. It is recommended to use rubber dam.

2. Prior to applying the filling material to the cavity the application of an adequate adhesive system is required (e.g. GLUMA® Comfort Bond, GLUMA® Solid Bond, iBOND® Self Etch by Heraeus Kulzer). Avoid contamination of saliva.

3. Class I–II and V-cavities: Apply CHARISMA® to the cavity in thin layers (max. 2 mm) and adapt carefully to the cavity walls. In case of considerable discolorations cavity linings with CHARISMA® opaque colours might become necessary.

4. Class III–IV-cavities: In case of anterior tooth cavities the setting-up of layers with an opaque nucleus often becomes necessary. CHARISMA® opaque Colours are especially suitable for this purpose.

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