**GC Fuji IX GP**RADIOPAQUE POSTERIOR GLASS IONOMER  
RESTORATIVE CEMENT

For use only by a dental professional in the recommended indications.

**RECOMMENDED INDICATIONS**

1. Class I and II restorations in deciduous teeth.
2. Non-load bearing Class I and Class II restorations in permanent teeth.
3. Intermediate restorative and base material for heavy stress situation in Class I and Class II cavities using sandwich laminate technique.
4. Class V and root surface restorations.
5. Core build-up.

**CONTRAINDICATIONS**


1. Pulp capping.
2. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

**DIRECTIONS FOR USE**

Powder / Liquid Ratio (g/g)	3.6 / 1.0
Mixing Time (sec.)	25-30"
Working Time (min., sec.)	2'00"
Net Setting Time (min., sec.)	2'20"
Test conditions : Temperature (23 +/- 1°C) Relative humidity (50 +/- 5%)	ISO 9917-1 : 2003 (E) (Glass polyalkenoate cement)

**1. POWDER AND LIQUID DISPENSING**

- a) The standard powder to liquid ratio is 3.6g/1.0g. (1 level scoop of powder to 1 drop of liquid).
- b) For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
- c) Hold the liquid bottle vertically and squeeze gently.
- d) Close bottles tightly immediately after use.

**2. MIXING**

Dispense powder and liquid onto the pad. Using the plastic spatula, divide the powder into 2 equal parts. Mix the first portion with all the liquid for 10 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-20 seconds (Fig. 1).

**3. RESTORATIVE TECHNIQUE**

- a) Prepare the tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping, use calcium hydroxide.
- b) Apply GC CAVITY CONDITIONER (10 seconds) or GC DENTIN CONDITIONER (20 seconds) to the bonding surfaces using a cotton pellet or sponge.
- c) Rinse thoroughly with water. Dry by blotting with a cotton pellet or gently blowing with an air syringe. DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening) (Fig. 2).
- d) Mix the required amount of GC Fuji IX GP. Working time is 2 minutes from the start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
- e) Transfer cement to the preparation using a syringe or other suitable instrument. Avoid air bubbles.
- f) Form the preliminary contour and cover with a matrix if required (Fig. 3).
- g) When set, immediately apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure) (Fig. 4).

**4. FINISHING**

- a) Final finishing under water spray using standard techniques can begin at 6 minutes after starting the mix (Fig. 5).
- b) Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.) (Fig. 6).
- c) Instruct the patient not to apply pressure for 1 hour.

**STORAGE**

Store at temperature of 4-25°C (39.2-77.0°F).

(Shelf life : 3 years from date of manufacture).

**SHADES**

A2, A3, A3.5, B2, B3, C4

Shade numbers according to Vita® shade guide.

Vita® is a trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

**PACKAGES**

- 1-1 package : 15g powder, 8g (6.4mL) liquid, powder scoop, mixing pad (No. 22), plastic spatula.
2. Bottle of 15g powder with scoop.
3. Bottle of 8g (6.4mL) liquid.

**CAUTION**

1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. DO NOT mix powder or liquid with any other glass ionomer components.

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Prior to use, carefully read  
the instructions for use. EN

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### RECOMMENDED INDICATIONS

- Class I and II restorations in deciduous teeth.
- Non-load bearing Class I and Class II restorations in permanent teeth.
- Intermediate restorative and base material for heavy stress situation in Class I and Class II cavities using sandwich laminate technique.
- Class V and root surface restorations.
- Core build-up.

### CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

### DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	3.6 / 1.0
Mixing Time (sec.)	25-30*
Working Time (min., sec.)	2'00"
Net Setting Time (min., sec.)	2'20"
Test conditions : Temperature (23 +/-1°C), Relative humidity (50 +/-5%)	

ISO 9917-1 : 2003 (E) (Glass polyalkenoaten cement)

#### 1. POWDER AND LIQUID DISPENSING

- The standard powder to liquid ratio is 3.6g/1.0g at level scoop of powder to 1.0ml of liquid.
- For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
- Hold the liquid bottle vertically and squeeze gently.
- Close bottles tightly immediately after use.

Fig.1

Fig.2

- Dispense powder and liquid onto the pad. Using the plastic spatula, divide the powder into 2 equal parts. Mix the first portion with all the liquid for 10 seconds to incorporate the remaining portion and mix the hole thoroughly for 15-20 seconds (Fig. 1).

#### 3. RESTORATIVE TECHNIQUE

- Prepare the tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping, use calcium hydroxide.
- Apply GC CAVITY CONDITIONER (10 seconds) or GC DENTIN CONDITIONER (20 seconds) to bonding surfaces using a cotton pellet or liquid blowing with an air syringe. DO NOT DESCALATE. Best results are obtained when prepared surfaces appear moist (glistening) (Fig. 2).
- Mix the required amount of GC Fuji IX GP. Working time is 2 minutes from the start of mixing at 23 °C (73.4°F). Higher temperatures will shorten working time.
- Transfer cement to the preparation using a syringe or other suitable instrument. Avoid air bubbles.
- Priming preliminary contour and cover with a matrix if required (Fig. 3).
- When set, immediately apply GC Fuji VARNISH (munk dry) or GC Fuji COAT LC (light cure) (Fig. 4).

#### 4. FINISHING

- Final finishing under water spray using standard techniques can begin at 6 minutes after starting the mix (Fig. 5).
- Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.)
- Instruct the patient not to apply pressure for 1 hour.

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\*Vita is a trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

#### PACKAGES

- 1-1 package : 15g powder, 8g (6.4mL) liquid, powder scoop, mixing pad (No. 22), plastic spatula.
2. Bottle of 8g with scoop.
3. Bottle of 8g (6.4mL) liquid.

#### CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- DO NOT mix powder or liquid with any other glass ionomer components.

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©Vita is een gedeponeerd handelsmerk van de Vita Zahnfabrik, Bad Sackingen, Duitsland
**GC GERMANY GmbH**
**GC AUSTRIA GmbH**