
## Directions For Use

# Identium®

Vinylsiloxanether® Impression Material ISO 4823


Gebrauchsinformation	Deutsch..... 2	Brugsanvisning	Svenska ..... 99
Directions for use	English ..... 17	Brugsinformation	Dansk..... 119
Mode d'emploi	Français ..... 31	Bruksanvisning	Norsk ..... 133
Istruzioni per l'uso	Italiano ..... 45	Kullanım Kılavuzu	Türkçe ..... 147
Instrucciones de uso	Español ..... 59	Instrukcja uży	Polski ..... 161
Instruções de uso	Português ..... 73	Gebruiksaanvatic	Nedderlands ..... 175
<b>Информация по использованию</b>	<b>Русский..... 87</b>		

---

## Directions For Use

# Identium®

Vinylsiloxanether® Impression Material ISO 4823

---

Made in Germany  
37402/0513

# Inhalt/Content/Contenu/Contenuto/Contenido/Conteúdo/содержимое/Innehåll

<b>Deutsch</b> .....	<b>3</b>	<b>Español</b> .....	<b>59</b>
Produktbeschreibung .....	4	Descripción del producto .....	60
Produktüberblick .....	7	Descripción general de productos .....	63
Technische Daten .....	8	Datos técnicos .....	64
Legende .....	16	Leyenda .....	72
<b>English</b> .....	<b>17</b>	<b>Português</b> .....	<b>73</b>
Product Description .....	18	Descrição do produto .....	74
Product Overview .....	21	Vista geral dos produtos .....	77
Technical Data .....	22	Dados técnicos .....	78
Legend .....	29	Legenda .....	86
<b>Français</b> .....	<b>31</b>	<b>Русский</b> .....	<b>87</b>
Description du produit .....	32	Описание продукта .....	88
Aperçu du produit .....	35	Обзор материалов .....	91
Caractéristiques techniques .....	36	Технические характеристики .....	92
Légende .....	44	Подписи к рисункам .....	98
<b>Italiano</b> .....	<b>45</b>	<b>Svenska</b> .....	<b>99</b>
Descrizione del prodotto .....	46	Produktbeskrivning .....	100
Panoramica dei prodotti .....	49	Produktinformation .....	103
Specifiche tecniche .....	50	Tekniska data .....	104
Legenda .....	58	Teckenförklaring .....	111

**Identium**<sup>®</sup> Heavy  
**Identium**<sup>®</sup> Heavy Fast  
**Identium**<sup>®</sup> Medium  
**Identium**<sup>®</sup> Medium Fast  
**Identium**<sup>®</sup> Medium Soft  
**Identium**<sup>®</sup> Light  
**Identium**<sup>®</sup> Light Fast  
**Identium**<sup>®</sup> Scan Heavy  
**Identium**<sup>®</sup> Scan Medium  
**Identium**<sup>®</sup> Scan Light

Hersteller:  
**Kettenbach GmbH & Co. KG**  
Im Heerfeld 7  
35713 Eschenburg, Germany  
[www.kettenbach.com](http://www.kettenbach.com)

Vertrieb in den USA:  
**Kettenbach LP**  
400 Oser Ave., Suite 1650  
Hauppauge, NY 11788, USA

**Vinylsiloxanether<sup>®</sup> Impression Material**

**ISO 4823**

Made in Germany  
37402/0513

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## Vinylsiloxanether<sup>®</sup> Impression Material ISO 4823

Made in Germany  
37402/0513

For professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.

## Product Description

Identium® products are addition-curing, elastomeric impression materials with parts of vinylpolysiloxane and polyether (Vinylsiloxanether®). They convince with excellent flow and hydrophilic properties, high tear strength, dimensional accuracy, and resistance to permanent deformation.

The Identium® product range consists of three different viscosities (heavy-bodied, medium-bodied, light-bodied), available in two delivery systems: standard 1:1 (50 ml automix cartridges) and 5:1 (362 ml foil bags) for use in most automatic dispensing and mixing systems. The Identium® products are all available in a Regular and Fast-set version.

Identium® Scan products (Medium, Heavy and Light) have been specifically developed for optical scanning using the CAD/CAM process (digitization; optimized for scanner systems working with infrared light) and guarantee the creation of impressions that offer the best scanning results (without using scanning powder). The applications are compatible with the impression techniques describe below.

## Indications/Techniques

**Identium® Heavy/Identium® Heavy Fast is to be used as a heavy-bodied material in one-step technique (double mix) for:**

- Impressions for crowns/bridges, inlays/onlays and veneer preparations
- Functional impressions
- Impressions for full or partial dentures
- Implant impressions

**Identium® Medium/Identium® Medium Fast is to be used as a medium-bodied tray or syringeable impression material in one-step technique (monophase or double mix) for:**

- Impressions for crowns/bridges, inlays/onlays and veneer preparations
- Implant impressions
- Fixation impressions
- Functional impressions

**Identium® Medium Soft is to be used as a medium-bodied tray or syringeable impression material in one-step technique (monophase or double mix) for:**

- Impressions for crowns/bridges, inlays/onlays and veneer preparations
- Functional impressions
- Reline impressions

**Identium® Light/Identium® Light Fast is to be used as syringeable impression material in one-step technique (double mix) for:**

- Impressions for crowns/bridges, inlays/onlays and veneer preparations
- Reline impressions
- Impressions for full or partial dentures

## Warnings

Do not use Identium® impression materials for two-step putty-wash impression technique or temporary reliner.

Do not use Identium® impression materials with other elastomeric impression materials.

Identium® Scan products (Medium, Heavy, Light) have been optimized for optical scanning during the CAD/CAM processes using scanner systems that work with infrared light. If other scanning technologies are used, optimum digital imaging cannot be guaranteed.

## Cautions

Do not use after expiration date.

Do not leave any residual impression material in the sulcus or oral cavity.

Do not swallow impression material! If swallowed: Consult a medical doctor if problems arise or persist.

Avoid contact with the eyes. If there is accidental contact with the eyes: Rinse immediately and thoroughly with an eye wash or water. Consult a medical doctor if problems arise or persist.

The products do normally not cause allergic reactions; however, for sensitive persons, an allergic reaction to the products cannot be ruled out.

Use of products containing active sulfur, aluminium chloride or nitrogen compounds (retraction cords containing ferric sulfate, polysulfide impression materials, etc.) in conjunction with this product will interfere with the setting reaction of the Vinylsiloxanether® material. Use of

these materials requires the area to be rinsed thoroughly to remove all residue. Avoid contact with latex gloves.

When taking impressions of areas with severe undercuts and wide interdental spaces, standard blocking-out techniques must be used.

When using a custom impression tray, ensure that sufficient space remains between the side of the tray and the teeth/jaw.

Do not allow the material to enter the sewer or water system to avoid environmental hazard.

Avoid contact with clothing, as the material cannot be removed by dry cleaning.

Do not use scanning powder when scanning impressions using Identium® Scan products (Medium, Heavy or Light).

## Product Overview

Product name	ISO 4823	Consistency (approx.) mm	Mixing ratio and total content volume per unit	Mixing technique	Mixing element
Identium® Heavy	Type 1, Heavy-bodied	32	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Identium® Heavy	Type 1, Heavy-bodied	34	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static green mixing tip MB Ø 6.5 mm
Identium® Heavy Fast	Type 1, Heavy-bodied	32	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Identium® Heavy Fast	Type 1, Heavy-bodied	34	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static green mixing tip MB Ø 6.5 mm
Identium® Medium	Type 2, Medium-bodied	36	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Identium® Medium	Type 2, Medium-bodied	36	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static green mixing tip MB Ø 6.5 mm
Identium® Medium Fast	Type 2, Medium-bodied	36	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Identium® Medium Fast	Type 2, Medium-bodied	36	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static green mixing tip MB Ø 6.5 mm
Identium® Medium Soft	Type 2, Medium-bodied	36	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Identium® Light	Type 3, Light-bodied	42	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static yellow mixing tip MB Ø 4.2 mm
Identium® Light Fast	Type 3, Light-bodied	42	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static yellow mixing tip MB Ø 4.2 mm

## Product Overview

Product name	ISO 4823	Consistency (approx.) mm	Mixing ratio and total content volume per unit	Mixing technique	Mixing element
Iidentium® Scan Heavy	Type 1, Heavy-bodied	32	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Iidentium® Scan Medium	Type 2, Medium-bodied	36	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Iidentium® Scan Light	Type 3, Light-bodied	42	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static yellow mixing tip MB Ø 4.2 mm

## Technical Data

Product name	Mixing ratio	Working time at 23° C/ 74° F ≤	Working time at 35° C/ 95° F ≤ <sup>1)</sup>	Intraoral setting time at 35° C/ 95° F ≥	Total setting time* ≥	Hardness (approx.) Shore	Linear dimensional change (approx.) %	Elastic-Recovery test (approx.) %	Strain in compression (approx.) %
Iidentium® Heavy	5:1	2 minutes	not applicable	3 minutes 30 seconds	5 minutes 30 seconds	A 60	-0.20	99.0	2.8
Iidentium® Heavy	1:1	2 minutes	not applicable	3 minutes 30 seconds	5 minutes 30 seconds	A 60	-0.20	99.0	2.5
Iidentium® Heavy Fast	5:1	1 minute 15 seconds	not applicable	2 minutes 15 seconds	3 minutes 30 seconds	A 60	-0.20	99.0	2.8
Iidentium® Heavy Fast	1:1	1 minute 15 seconds	not applicable	2 minutes 15 seconds	3 minutes 30 seconds	A 60	-0.20	99.0	2.5
Iidentium® Medium	5:1	2 minutes	1 minute 20 seconds	3 minutes 30 seconds	5 minutes 30 seconds	A 60	-0.20	99.0	2.3
Iidentium® Medium	1:1	2 minutes	1 minute 20 seconds	3 minutes 30 seconds	5 minutes 30 seconds	A 60	-0.20	99.0	2.3

## Technical Data

Product name	Mixing ratio	Working time at 23° C/ 74° F ≤	Working time at 35° C/ 95° F ≤ <sup>1)</sup>	Intraoral setting time at 35° C/ 95° F ≥	Total setting time* ≥	Hardness (approx.) Shore	Linear dimensional change (approx.) %	Elastic-Recovery test (approx.) %	Strain in compression (approx.) %
Identium® Medium Fast	5:1	1 minute 15 seconds	40 seconds	2 minutes 15 seconds	3 minutes 30 seconds	A 60	-0.20	99.0	2.3
Identium® Medium Fast	1:1	1 minute 15 seconds	40 seconds	2 minutes 15 seconds	3 minutes 30 seconds	A 60	-0.20	99.0	2.3
Identium® Medium Soft	5:1	2 minutes	1 minute 20 seconds	3 minutes 30 seconds	5 minutes 30 seconds	A 50	-0.20	99.0	3.3
Identium® Light	1:1	2 minutes	1 minute 20 seconds	3 minutes 30 seconds	5 minutes 30 seconds	A 46	-0.20	99.0	3.8
Identium® Light Fast	1:1	1 minute 15 seconds	40 seconds	2 minutes 15 seconds	3 minutes 30 seconds	A 46	-0.20	99.0	3.8
Identium® Scan Heavy	5:1	2 minutes	not applicable	3 minutes 30 seconds	5 minutes 30 seconds	A 60	-0.20	99.0	2.8
Identium® Scan Medium	5:1	2 minutes	50 seconds	3 minutes 30 seconds	5 minutes 30 seconds	A 60	-0.20	99.0	2.3
Identium® Scan Light	1:1	2 minutes	1 minute 20 seconds	3 minutes 30 seconds	5 minutes 30 seconds	A 46	-0.20	99.0	3.8

\*Total setting time (removal time from mouth) from start of mix.

<sup>1)</sup>With a reduced working time the total setting time of 3 minutes 30 seconds for the Fast set version and 5 minutes 30 seconds for the regular set version has always to be heeded.

## Note

To ensure optimal impressions the temperature of the material should not deviate from 23° C (74° F) before applying. Otherwise, working and setting times will be affected.

## Impression tray: Preparation and adhesives

In principle, all common impression trays can be used if a relevant dynamic pressure is guaranteed. When retention is not sufficient, to ensure a strong bond to the impression material, brush the impression tray with a thin film of Identium® Adhesive prior to loading the tray with impression material. Allow to dry as per instructions for use and cover with a sufficient amount of impression material.

## Preparation of the Plug & Press® foil bag system

### Foil bag material for use in most automatic dispensing and mixing systems

Before using the base foil for the first time, remove the white safety pin from the activating head of the large foil bag by turning it in the direction of the arrow and pulling it out (Figure 1, Page 112/113).

Place the large foil bag with the integrated activating head in the cartridge body. Ensure that the notches on the activating head and cartridge case are aligned (Figure 2, Page 112/113).

Firmly press the activating head into the final position on the cartridge body. The foil bag is automatically pierced by a pin when the activating head is pressed down (Figure 3, Page 112/113).

Follow the same procedure as in Figure 2 and 3 for the smaller foil bag. (Note: There is no safety pin on the catalyst foil bag.)

The use of dynamic mixers and cartridge bodies supplied by Kettenbach is a prerequisite for an optimum result.

## Application using the Plug & Press® Dispenser in combination with the Kettenbach's dynamic mixing tip

- Move the plungers to the top by turning the control knob (Figure 4, Page 112/113).
- Insert the cartridge body into the unit (Figure 5, Page 112/113).
- Turn the control knob to move the plungers into the cartridge body and continue turning until the plungers

come into contact with the foil bags (Figure 6, Page 112/113).

- Once there is contact, press one of the two start keys (at the front on top of the unit) to dispense the material (Figure 7, Page 112/113).
- Before fitting the dynamic mixer, dispense material until equal amounts of base and catalyst appear (Figure 7, Page 112/113), wipe excess. Firmly place dynamic mixer (Figure 8, Page 112/113) and lower locking lever.

Load the tray with the required amount of material. Hold the impression tray at an angle and press lightly against the tray (Figure 9, Page 112/113). Leave the dynamic mixer in the material when dispensing. Leave the filled dynamic mixer on the cartridge body as a seal.

Before next use, release the locking clip to remove the used dynamic mixer and check that the outlets in the activating heads are not blocked. Fit a new dynamic mixer, lower locking lever and continue as usual.

When finished simply remove the empty foil bags and activating heads from the cartridge body and discard them. The cartridge body is reusable (Figure 10, Page 112/113).

The activating heads are disposable.

If not using a Kettenbach Plug & Press® Dispenser, please follow the manufacturer's instructions for your automatic dispensing and mixing system.

## **Application using the Sympress Dispenser in combination with the Kettenbach's dynamic mixing tip**

- Allow the plungers (with cover closed) to return to the original position by pressing the plunger return button (Figure 1, Page 114/115). Next, open the cover (Figure 2, Page 114/115).
- Insert the cartridge body with the foil bags into the device (if foil bag material is already in use, the filled mixer will be on the cartridge body as the closure) (Figures 2 and 3, Page 114/115).
- Mount the dynamic mixer in place (if you are already using foil bag material, first remove the mixer, which is being used as the closure and replace it with the new dynamic mixing tip) (Figure 4, Page 114/115).
- Close the locking lever (Figure 5, Page 114/115).

- Close the device cover (Figure 6, Page 114/115).
- Dispense the material, dispose of the initial approx. 3 cm (until a homogenous mix of base and catalyst appear; see coloring of the final impression material) (Figures 7 and 8, Page 114/115).
- Load the tray (or application syringe) (Figure 9, Page 114/115). As you do, hold the impression tray at an angle and press lightly against it. Leave the mixing tip in the material while dispensing to prevent air from getting trapped.
- Leave the filled dynamic mixer in place as a seal.

Before next use, release the locking clip to remove the used dynamic mixer and check that the outlets in the activating heads are not blocked. Fit a new dynamic mixer, lower locking lever and continue as usual.

When finished simply remove the empty foil bags and activating heads from the cartridge body and discard them. The cartridge body is reusable (Figure 10, Page 114/115). The activating heads are disposable.

If you are using an automated dosage and mixing system other than the Sympress Dispenser, please follow the instructions for your respective system.

## **Instructions for use:**

### **Cartridge material for use in Applyfix® 4 dispensing gun DS-50 1:1/2:1**

Insert the cartridge into the Applyfix® 4 dispensing gun DS-50 1:1/2:1. Ensure that the notches on the cartridge base are pointing downwards. The clasp will not close if the cartridge is not properly inserted into the dispensing gun.

Remove the cartridge cap. The cap can be replaced after initial use.

Advance plunger into cartridge and dispense a small amount of impression material until equal amounts are extruded at the same rate.

Install a mixing tip on the cartridge and turn cap 90 degrees in a clockwise direction to lock in place.

Inject the required amount of material, e.g., directly into the tray or the preliminary impression, or into a delivery syringe. To apply the material around the prepared teeth use the delivery syringe or the cartridge with an intraoral tip on the mixer. Check to ensure materials are set prior to mouth removal.

Leave used mixing tip on the cartridge after use or replace

the cartridge cap. Before using the cartridge again, remove cartridge cap or remove and discard the used mixing tip.

Check the cartridge openings to ensure no polymerized material is present. Should this be the case, unblock the openings and dispense a small amount of impression material until equal amounts are extruded at the same rate.

Install a new mixing tip and continue as described above.

The Applyfix® 4 dispensing gun DS-50 1:1/2:1 can be sterilized in the autoclave.

## **Instructions for use:** **Applyfix® 5 application syringe**

Insert a new syringe tip into tip retainer. Screw tip retainer with the syringe tip onto the syringe cylinder. (Figure 2, Page 116)

If necessary, lubricate the ring washer very slightly with petrolatum for an easy action of the syringe piston.

Prepare the preferred impression material (the Applyfix® 5 application syringe is indicated for light-bodied or medium-bodied consistency impression material only).

Attach the open syringe cylinder of the Applyfix® 5 directly to the mixing tip of the preferred foil bag impression material. (Figure 3, Page 116)

Fill the needed amount of material directly into the syringe cylinder. Insert the syringe piston into the syringe cylinder (Figure 4, Page 116) and apply the impression material into the patient's mouth as desired.

### Cleaning

Clean the Applyfix® 5 application syringe using with the brush (Figure 5, Page 116) after using and removing of the full hardened impression material. The Applyfix® 5 can be reused and autoclaved several times and should be replaced once it shows signs of mechanical damage. Do not reuse the syringe tips, to avoid contamination.

### Cautions

Only use one material per Applyfix® 5 application syringe. Contamination of different impression materials renders the impression material unusable.

### **Note**

Keep in the original container.  
Selected Kettenbach materials/products are available in certain markets only.

## **Storage conditions:**

Store in a dry place at room temperature; do not expose to sunlight and avoid exposure to extreme temperature fluctuations.

## **Disinfection:**

Impressions can be disinfected using for example, a 2% acidic glutaraldehyde solution. Use disinfection solutions-specific for impression materials (e.g., Silosept®). See manufacturer instructions.

## **Electroplating:**

Impressions may be silver-plated with an alkaline-plating solution only.

## **Model preparation:**

Impressions with Identium® may be poured immediately after disinfection of the impression up to several weeks with standard dental stone (class III - V). Impressions can also be poured several times.

## **Trademarks**

Applyfix® 4, Applyfix® 5, Identium®, Silosept® and Plug & Press® Dispenser are registered trademarks of Kettenbach GmbH & Co. KG.

## **Limitation of liability**

Except where prohibited by law, Kettenbach GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.


The information provided for Kettenbach products is based on comprehensive research and experience in application technology.

Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development.

However, users must comply with and consider all recommendations and information in connection with any use.

Updated: 28 January 2013

## Legend

Conformity with MDD 93/42/EEC		Use by	
Temperature limitation		Keep dry	
Keep away from sunlight		Catalog number	<b>REF</b>
Batch code		Caution, consult accompanying documents	
Dual system (Germany only)		Diameter	
Cartridge		Static mixing tip MB (single use)	
Mixing tip Type B	<b>MB</b>	Intraoral tip (single use)	
Adhesive		Do not resterilize	
Cartridge body		Foil bags	
Dynamic mixer (single use)		Rx only For professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.	
Applyfix® 5 application syringe		Applyfix® 5 syringe tips (single use)	
Millimeter	<b>mm</b>	Milliliter	<b>ml</b>
Less than or equal to	<b>≤</b>	Greater than or equal to	<b>≥</b>

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