
Silginat®

Silginat® Strawberry

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**Vinyl polysiloxane impression material as an
alternative to alginate applications**

ISO 4823

Made in Germany
37406/3412

For professional use only. Caution: Federal (U.S.A.) Law
restricts this device to sale by or on the order of a dentist.

Vinyl Polysiloxane Impression Material ISO 4823

Product Description

Alginate Alternative: Silginat® and Silginat® Strawberry are addition-curing, elastomeric impression materials intended to be used as an alternative to traditional alginate materials (e.g., for anatomical impressions). Available in 1:1 38 ml cartridges (optimal volume for a single arch or two quadrant impressions) suitable for the Kettenbach Applyfix® 4 dispensing gun and in 5:1 362 ml Plug and Press® foil bags for use in most commercially available mixing machines, e.g., Sympress.

Indications

Silginat® and Silginat® Strawberry are useful as alginate replacements for making impressions for:

- Opposing dentitions
- Fabricating temporary crowns and bridges
- Anatomical models
- Fabricating simple removable prosthetic restorations
- Producing models for the construction of trays for whitening solutions, bites and surgical guides.

- Orthodontic work
- Case study models
- Fabricating removable retainers and splints

Cautions

Do not use after expiration date.

Do not leave any residual impression material in the sulcus or oral cavity.

Do not swallow impression material! If swallowed: consult a medical doctor if problems arise or persist.

Avoid contact with eyes. If accidental contact with the eyes occurs, rinse immediately and thoroughly with an eye wash or water. Consult a medical doctor if problems arise or persist.

The product does not normally cause allergic reactions. However, for sensitive persons, an allergic reaction to the product cannot be ruled out.

Use of products containing active sulfur, aluminium chloride or nitrogen compounds (retraction cords containing ferric sulfate, polysulfide impression materials, etc.) in conjunction with this product will interfere with the setting reaction of the vinyl polysiloxane material. Use of the materials requires the area to be rinsed thoroughly to remove all residue, before an impression is taken. Avoid contact with latex gloves.

Do not allow the material to enter the sewer or water system to avoid environmental contamination.

Avoid contact with clothing, since the material cannot be removed by dry cleaning.

Product Overview

Product name	ISO 4823	Consistency (approx.) mm	Mixing ratio and total content volume per unit	Mixing technique	Mixing element
Silginat®	Type 2, Medium-bodied	33	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Silginat®	Type 2, Medium-bodied	35	1:1, 38 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static green mixing tip MB Ø 6.5 mm
Silginat® Strawberry	Type 2, Medium-bodied	33	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Silginat® Strawberry	Type 2, Medium-bodied	35	1:1, 38 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static green mixing tip MB Ø 6.5 mm

Technical Data

Product name	Mixing ratio	Working time at 23° C/74° F ≤	Intraoral setting time at 35° C/95° F ≥	Total setting time* ≥	Hardness (approx.) Shore	Linear dimensional change (approx.) %	Elastic-Recovery test (approx.) %	Strain in compression (approx.) %
Silginat®	5:1	1 minute 30 seconds	1 minute 30 seconds	3 minutes	A 42	- 0.20	99.5	7.2
Silginat®	1:1	1 minute 15 seconds	1 minute 15 seconds	2 minutes 30 seconds	A 42	- 0.20	99.5	7.2
Silginat® Strawberry	5:1	1 minute 30 seconds	1 minute 30 seconds	3 minutes	A 42	- 0.20	99.5	7.2
Silginat® Strawberry	1:1	1 minute 15 seconds	1 minute 15 seconds	2 minutes 30 seconds	A 42	- 0.20	99.5	7.2

*Total setting time (removal time from mouth) from start of mix.

Note

To ensure optimal impressions, the temperature of the material should not deviate from 23° C (74° F) before applying. Otherwise, working and setting times will be affected.

Impression tray: Preparation and adhesives

In principle, all common impression trays can be used if a relevant dynamic pressure is guaranteed. When retention is not sufficient, ensure a strong bond to the impression material, by brushing the impression tray with a thin film of adhesive for addition-curing silicone prior to loading the tray with impression material (e.g. Panasil® Adhesive). Allow to dry per manufacturer's instructions.

Preparation of the foil bag system

Foil bag material for use in automatic dispensing and mixing systems

Before using the base material for the first time, remove the white safety pin from the activating head of the large foil bag by turning it in the direction of the arrow and pulling it out (Figure 1, Page 74/75).

Place the large foil bag with the integrated activating head

in the cartridge body. Ensure that the notches on the activating head and cartridge case are aligned (Figure 2, Page 74/75).

Firmly press the activating head into the final position on the cartridge body. The foil bag is automatically pierced by a pin when the activating head is pressed down (Figure 3, Page 74/75).

Follow the same procedure as in Figure 2 and 3 for the catalyst material. (Note: There is no safety pin on the catalyst foil bag.)

The use of dynamic mixers and cartridge bodies supplied by Kettenbach is a prerequisite for an optimum result.

Application using the Plug & Press® Dispenser in combination with the Kettenbach's dynamic mixing tip

- Move the plungers to the top by turning the control knob (Figure 4, Page 74/75).
- Insert the cartridge body into the unit (Figure 5, Page 74/75).
- Turn the control knob to move the plungers into the cartridge body, and continue turning until the plungers come into contact with the foil bags (Figure 6, Page 74/75).
- Once there is contact, press one of the two start keys (at the front on top of the unit) to dispense the material (Figure 7, Page 74/75).
- Before fitting the dynamic mixer, dispense material until equal amounts of base and catalyst appear (Figure 7, Page 74/75); wipe excess. Firmly place dynamic mixer (Figure 8, Page 74/75) and lower locking lever.

Load the tray with the required amount of material. Hold the impression tray at an angle and press lightly against

the tray. Leave the dynamic mixer in the material while dispensing (Figure 9, Page 74/75). Leave the filled dynamic mixer on the cartridge body as a seal.

Before next use, release the locking clip to remove the used dynamic mixer and check that the outlets in the activating heads are not blocked. Fit a new dynamic mixer, lower locking lever and continue as usual.

When finished, simply remove the empty foil bags and activating heads from the cartridge body and discard them. The cartridge body is reusable (Figure 10, Page 74/75). The activating heads are disposable.

If not using a Kettenbach Plug & Press® Dispenser, please follow the manufacturer's instructions for your automatic dispensing and mixing system.

Application using the Sympress Dispenser in combination with the Kettenbach's dynamic mixers

- Allow the plungers (with cover closed) to return to the original position by pressing the plunger return button (Figure 1, Page 76/77). Next, open the cover (Figure 2, Page 76/77).

- Insert the cartridge body with the foil bags into the device (if foil bag material is already in use, the filled mixer will be on the cartridge body as the closure) (Figures 2 and 3, Page 76/77).
- Mount the dynamic mixer in place (if you are already using foil bag material, first remove the mixer, which is being used as the closure and replace it with the new dynamic mixer). (Figure 4, Page 76/77).
- Close the locking lever (Figure 5, Page 76/77).
- Close the device cover (Figure 6, Page 76/77).
- Dispense the material, dispose of the initial approx. 3 cm (until a homogenous mix of base and catalyst appear; see coloring of the final impression material) (Figures 7 and 8, Page 76/77).
- Load the tray (or application syringe) (Figure 9, Page 76/77). As you do, hold the impression tray at an angle and press lightly against it. Leave the mixing tip in the material while dispensing to prevent air from getting trapped.

- Leave the filled dynamic mixer in place as a seal.

Before next use, release the locking clip to remove the used dynamic mixer and check that the outlets in the activating heads are not blocked. Fit a new dynamic mixer, lower locking lever and continue as usual.

When finished simply remove the empty foil bags and activating heads from the cartridge body and discard them. The cartridge body is reusable (Figure 10, Page 76/77). The activating heads are disposable.

If you are using an automated dosage and mixing system other than the Sympress Dispenser, please follow the instructions for your respective system.

Instructions for use: Cartridge material for use in Applyfix® 4 dispensing gun DS-50 1:1/2:1 (Page 78)

Insert the cartridge into the Applyfix® 4 dispensing gun DS-50 1:1/2:1. Ensure that the notches on the cartridge base are pointing down. The clasp will not close if the cartridge is not properly inserted into the dispensing gun.

Remove the cartridge cap. The cap can be replaced after initial use.

Advance plunger into cartridge and dispense a small amount of impression material until equal amounts are extruded at the same rate.

Install a mixing tip on the cartridge, and turn cap 90 degrees clockwise to lock in place.

Inject the required amount of material, directly into the tray.

Leave used mixing tip on the cartridge after use or replace the cartridge cap. Before using the cartridge again, remove cartridge cap or remove and discard the used mixing tip.

Check the cartridge openings to ensure no polymerized material is present. Should this be the case, unblock the openings and dispense a small amount of impression material until equal amounts are dispensed at the same rate.

Install a new mixing tip and continue as described above.

The Applyfix® 4 dispensing gun DS-50 1:1/2:1 can be sterilized in the autoclave.

Storage conditions:

Store in a dry place at room temperature; do not expose to sunlight and avoid exposure to extreme temperature fluctuations.

Disinfection:

Impressions can be disinfected using for example, a 2% acidic glutaraldehyde solution. Use disinfection solutions specific for impression materials (e.g., Silosept®). See manufacturer instructions.

Electroplating:

Impressions may be silver-plated with an alkaline-plating solution only.

Model preparation:

Impressions with Silginat® and Silginat® Strawberry may be poured immediately after disinfection of the impression up to several weeks with standard dental stone (e.g., Tewerock®, Tewestone®). Impressions can also be poured several times. A surfactant solution is not required.

Trademarks

Applyfix® 4, Silginat®, Silosept®, Panasil®, Tewaterock®, Tewestone® and Plug & Press® Dispenser are registered trademarks of Kettenbach GmbH & Co. KG.

Note

Keep in the original container.

Selected Kettenbach materials are available in certain markets only.

Limitation of liability

Except where prohibited by law, Kettenbach GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.


The information provided for Kettenbach products is based on comprehensive research and experience in application technology.

Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development.

However, users must comply with and consider all recommendations and information in connection with any use.

Updated: 30th January 2012

Legend

Conformity with MDD 93/42/EEC		Use by	
Temperature limitation		Keep dry	
Keep away from sunlight		Catalog number	REF
Batch code		Caution, consult accompanying documents	
Dual system (Germany only)		Diameter	
Cartridge		Static mixing tip MB (single use)	
Mixing tip Type B	MB		
Adhesive		Single use	
Cartridge body		Foil bags	
Dynamic mixer (single use)		Rx only For professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.	
Millimeter	mm	Milliliter	ml
Less than or equal to	≤	Greater than or equal to	≥

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