## D Gebrauchsinformation

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1 Kartusche à 76 g Paste  
15 Automix-Tips

## GB Instruction for use

DC Crown&Bridge Material is a self-curing composite for the fabrication of temporary crowns and bridges, inlays, onlays and veneers. The 2 component multifunctional methacrylic material is automatically mixed. The new fluorescence allows for naturally aesthetic, life-like temporaries. The low setting temperature (maximum 38° C/ 100° F), ensures that the sensible pulpa tissue is treated with care and postoperative sensibilities are reduced.

### Indication:

Temporary crowns, bridges, inlays, onlays, partial crowns, veneers and long-term temporaries

### Timing:

Insertion in the mouth:  
0:00 – 0:45 min  
Removal from the mouth:  
2:00 – 3:00 min  
End of setting:  
6:00 – 7:00 min

When using a cartridge for the first time, dispense a pea size amount of material and discard.

Leave the used cannula on the cartridge until next use! It serves as a cap.

### Recommended procedure:

- Prior to preparation or extraction of a tooth, take a quadart impression using alginate, silicone or a thermoplastic impression material.
  - Dry prepared teeth and lightly lubricate prepared teeth and surrounding tissue and eventual resin based core built-ups with Vaseline or a similar separating medium.
  - Match shade and place appropriate DC Crown&Bridge Material into the dispenser.
  - Dispense DC Crown&Bridge Material into the impression or matrix. Dispense the material first onto the occlusal surface, and then bring it gingivally. Make sure to overfill it slightly. To prevent bubbles, it is important to keep the tip of the mixing cannula immersed in the material.
  - Insert the impression or matrix, filled with DC Crown&Bridge Material within 45 seconds with gentle pressure over the prepared teeth and hold firmly in place.
  - After approximately 2-3 minutes, the provisional crown or bridge can be removed together with the impression from the prepared teeth.
- The setting reaction has to be monitored intraorally (e.g. with a scaler), as the temperature in the mouth has a significant influence on the setting time and removal of the provisional is only possible during its rubbery phase.**
- The crown or bridge may be contoured and polished after setting 6 additional minutes. We recommend removing the removal of the soft, sticky inhibition layer on the surface of the DC Crown&Bridge Material restoration (e.g. with DC Clean or alcohol). The temporary can be shaped with slow-speed acrylic burs and disks. Use appropriate safety measures. Use high-speed finishing diamond burs in the mouth to adjust occlusion and further correct embrasures.
  - Glaze surface with a light-curing high gloss varnish or polish manually with a rubber polisher or a goat's hair wheel.
  - Cement the temporary with any temporary cement (e.g. DC EugenolfreeCem).
- Cements containing eugenol may inhibit**

**polymerization of resin based luting cements for permanent cementation.**

### Please note:

- Carving in molar areas with absent teeth, be sure that the interproxial areas should create a ridge for pontic contact between the prepared teeth.
- The model is made from a negative image of an impression using a plastic template. This can be done using a vacuum formmachine or silicone compression technique to fabricate a plastic matrix.
- To further improve the precision of larger temporaries, the restoration can be reinserted into the impression until reaching its final hardness.

### Repair of temporaries:

Repair of an older temporary: Using a bur or air abrasion unit, slightly roughen the damaged surface. Apply some bonding agent according to manufacturer's instructions. For the repair of a freshly made provisional, a preparation of the surface or the use of a bonding agent is not necessary.

- Dispense a fresh mix of material onto all surfaces.
- Hold parts together for 3 minutes.
- Wait 6 minutes before removing excess with rotary instruments.

To speed up polymerization, the repaired temporary may be put into water of 50° C/122° F for a few minutes. The DC Crown&Bridge Material temporary can alternatively be repaired with the flowable lightcuring composite DC FlowFill. For information about the handling of the product, please refer to the instructions.

### Technical data:

Compressive strenght	200 MPa
Transverse strenght	60 MPa
Diametral tensile strenght	31 MPa
Barcol hardness	20
Water sorption	25 µg/mm <sup>3</sup>
Temperature rise during polymerization	~ 38° C/100° F

### Contact with skin and eyes should be avoided!

### Composition:

Glass filler in a matrix of multifunctional methacrylates; catalysts, stabilizers, additives. Free of methyl methacrylate and peroxides!

### Storage:

Store dry at maximal 25° C/77° F. Only use at room temperature (Refrigerated material is more viscous and cures more slowly). Do not use after expiry date.

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