WELLSACLEAN S+

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REF WCC106

REF WCC116

Produkttyp

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REF WCC116

Product description

WELLSACLEAN S+ SCHNELLDESINFEKTION is a ready to use, alcoholic solution for surface disinfection. The product must not be diluted by the user. It contains quarternary ammonium compounds and is free of aldehydes and phenols. After the contact time has elapsed, treated surfaces dry quickly. WELLSACLEAN S+ SCHNELLDESINFEKTION can be applied for wipe as well as spray disinfection.

Areas of application

WELLSACLEAN S+ SCHNELLDESINFEKTION is suitable for rapid disinfection of alcohol-resistant surfaces of non-invasive, non-critical medical devices, medical equipment as well as all kinds of wipeable surfaces.

Suitable for:

- Surfaces close to patients, dental chairs, examination couches, operating tables
- Areas in medical and dental practices which are used for medical purposes
- Rapid disinfection between surgeries and procedures
- Clothes like OR shoes, OR aprons, X-ray aprons
- Alcohol-resistant surfaces of medical equipment
- Materials according to material compatibility (see page 2)

Not tested for food and veterinary areas. Not suitable for disinfection of skin and hands. Not suitable for disinfection of invasive medical devices. Do not use for disinfection of semi-critical and critical medical devices.

Further information on cleaning and disinfection of surfaces and medical devices:

Recommendation of the Commission for Hospital Hygiene and Infection Prevention at the Robert-Koch Institute (RKI) entitled „Hygiene Requirements for Cleaning and Disinfecting Surfaces“. Available online.

Personal protection

Put on protective gloves, face mask and safety glasses prior to each disinfection. Provide adequate ventilation. Keep away from unprotected individuals. Avoid contact with eyes. Avoid breathing vapours. Keep away from sources of ignition.

Application

In case the surfaces are contaminated with organic material (e.g. blood, secretions etc.): Remove the visible contamination with a disposable cloth soaked with disinfectant and discard it. Subsequently, disinfect surface. Please note that the disinfection does not substitute the cleaning of the surface.

Wipe disinfection: Apply the ready to use solution to a clean, dry single use cloth or directly onto the surface. Wipe the surface applying mild pressure (wet wiping) and keep the surface completely moist during contact time.

Spray disinfection: Spray disinfection may only be applied for inaccessible, twisted surfaces. For spray disinfection, apply the product undiluted directly onto the surface. Ensure that the surface is completely moist and keep moist during contact time.

After the contact time has elapsed: Remove residues with a clean, dry, undyed disposable cloth or allow to dry.

Special advice

Read safety datasheet and technical datasheet before use.

Do not dilute the product. Do not mix with additives or other products.

For professional use only.

Do not use after expiry date.

Electronic devices must be disconnected from their power supply before disinfection. Ensure that the disinfectant does not seep into electronic devices.

WELLSACLEAN S+

SCHNELLDESINFEKTION

REF WCC106

REF WCC116

Efficacy spectrum

Efficacy (ready to use product, at room temperature)	Standard	Contact Time
BACTERIA		
S. aureus, E. hirae, P. aeruginosa, E. coli, MRSA → (dirty conditions)	EN 13727	30 sec
S. aureus, E. hirae, P. aeruginosa, E. coli → (dirty conditions*)	EN 13697	5 min
S. aureus, E. hirae, P. aeruginosa	EN 16615	1 min
S. aureus, E. hirae, P. aeruginosa → (dirty conditions)	Standard methods VAH*2	5 min
MYCOBACTERIA/TB		
M. terrae, M. avium → (dirty conditions)	EN 14348	1 min
YEASTS		
C. albicans → (dirty conditions)	EN 13624	30 sec
C. albicans → (dirty conditions)	EN 16615	1 min
C. albicans → (dirty conditions*)	EN 13697	5 min
C. albicans → (dirty conditions)	Standard methods VAH*2	5 min
VIRUSES		
Enveloped viruses		
MVA Modified Vacciniavirus Ankara → (dirty conditions)	EN 14476	30 sec
BVDV → (dirty conditions)		
Non-enveloped viruses		
Adenovirus type 5 → (dirty conditions)	EN 14476	5 min
Human Rotavirus → (clean conditions)		
MNV murine Norovirus → (dirty conditions)		
MNV murine Norovirus → (dirty conditions)	EN 16777	5 min
Adenovirus type 5 → (dirty conditions)		

*dirty conditions according to EN 13697:2015-06 Quantitative non-porous surface test for the evaluation of bactericidal and/or fungicidal activity of chemical disinfectants used in food, industrial, domestic and institutional areas.

*2 standard methods of the VAH (Association of Applied Hygiene) for surface disinfection by wiping (quantitative carrier test with mechanical action, practical conditions)

Material compatibility

Metals	aluminium, stainless steel, copper, brass
Plastics	PSU, PVC, silicone
Rubber	EPDM, natural rubber
Advice	Upon compliance to the intended use and the indicated contact time, no material damage is to be expected with the above mentioned materials. Avoid contact longer than 48 hours with linoleum flooring. Avoid contact with acrylates. May cause stress corrosion on acrylic glass. Test material compatibility for polyacrylates, polycarbonates as well as not mentioned materials in an inconspicuous place before use.

Composition

100 g ready to use solution contain:

33.66 g ethanol, 15.84 g 2-propanol, 0.05 g didecylmethyl-ammonium chloride 70%, perfume*³ (not applicable for perfume-free)

Physico-chemical data

appearance: clear and colourless
 pH-value: 5.0 – 6.0
 density (20 °C): 0.900 – 0.920
 flash point (EN 22719): 22.75 + 0.5 °C

Storage conditions

Keep container tightly closed and in a well-ventilated place. Keep cool. Protect from sunlight. Keep away from sources of ignition. No smoking. Take precautionary measures against static discharge. Store in original container. Store in an upright position. See safety datasheet for further advice.

Labelling according to Regulation (EC) No 1272/2008 (CLP)

For information on labeling according to Regulation (EC) No 1272/2008 (CLP) see safety datasheet and label.

UN 1987

Identification, certificates

Medical device class II a / MD 0108
 Certificate of efficacy issued by VAH
 (Association of Applied Hygiene)

CE 0123

Expert opinions and further information

WELLSACLEAN S+ SCHNELLDESINFEKTION has been assessed by independent laboratories.

To exclude adverse effects on the skin upon unintentional contact with the product, a dermatologically survey was carried out by the company Dermatest.

For further details concerning expert opinions, information on products and answers to specific questions contact us by e-mail (info@wellsamed.de).

The manufacturer has a manufacturing authorisation according as per §13 German medicines act as well as a GMP- and an EN 13485 certificate.

Please report serious incidents related to this product to the manufacturer and the respective authorities.

*3 Please note that perfume may cause allergies.

WELLSACLEAN S+

SCHNELLDESINFEKTION

REF WCC106

REF WCC116

Related medical devices manufactured with WELLSACLEAN S+ SCHNELLDESINFEKTION:

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(ready to use, alcohol-soaked disinfection wipes for application with dispenser bucket)

REF **WCW407** 70 sheets of fleece towel roll

REF **WCW409** 90 sheets of fleece towel roll

PLIWA Hygiene GmbH,
Pliwastraße 2, 34323 Malsfeld, Germany

Distributed by:
wellsamed GmbH
Brahestraße 16, 04347 Leipzig, Germany