# CHARISMA® *flow*

66034065/02 – 07/09

## CHARISMA® *flow*

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## CHARISMA® *flow*

Charisma flow is a flowable, light curing, radiopaque hybrid composite used for adhesive, tooth-coloured anterior and posterior restorations. Due to its low viscosity Charisma flow provides excellent coverage of the tooth structure.

**Composition:** Charisma Flow is based on multifunctional methacrylate monomers (EBADMA/TEGDMA) and contains approximately 62% by weight or 38% by volume inorganic fillers such as Ba-Al-F silicate glass and SiO<sub>2</sub>. The particle size of these fillers ranges from 0.005µm and 5µm.

**Shades:** A1, A2, A3, A3.5, B2, B3, OA2, OA4, OB2, SB, Baseline

**Applications:**

- Fissure sealing
- Distributed fissure sealing
- Class V fillings
- Minimally invasive Class I fillings
- Minimally invasive Class II fillings
- Minimally invasive Class III fillings
- Smaller contour and shade adjustment on the enamel and dentine
- Small surface repairs to direct and indirect restorations in combination with a suitable bonding agent
- Springing of teeth
- Cavity lining
- Bracket retention
- Bracket bonding

**Kind of application:**

- Prepare the tooth structure for treatment by cleaning it with a fluoride free polishing paste. Before treatment, while the teeth are still moist, select the Charisma flow shade using the Charisma shade guide.
- 1. Remove all caries and completely minimize invasive preparation of the tooth structure for adhesive fillings. Then clean and dry the cavity. The use of a rubber dam is recommended.
- 2. A quality adhesive system (e.g. Bond® Self Etch, GLUMA® Comfort® Bond + Desensitizer, GLUMA® Solid Bond) should be applied to the enamel and dentin in accordance with the manufacturer's instructions before insert- ing the filling material. Avoid any contamination of conditioned or adhesive-treated surfaces with moisture, blood or saliva.
- 3. Apply Charisma flow in the layers (max. 2mm, Baseline max. 1mm) to the cavity. When applying Charisma flow ensure that it flows slowly and evenly to prevent air bubbles – apply a gentle stream of air to facilitate maximum marginal contact. If necessary, remove any air bubbles with a probe before polymerisation.
- 4. Polymerise each layer according to the table below using a suitable light-curing unit (e.g. Translux® Power Blue®) for 20 or 40sec. A dispersion layer forms during polymerisation. This layer should not be touched or removed if other composite layers are to be applied.

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