FR

COMPOSITION

Chlorure d'aluminium hexahydraté 25 % m/V, oxyquinol, excipient hydroalcooolique.

PROPRIÉTÉS

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UK

COMPOSITION

Hexahydrated aluminium chloride 25 % w/V, oxyquinol, hydroalcoholic excipient.

PROPERTIES

Aluminium chloride has an astringent effect. It induces local haemostasis and the contraction of the upper strata of the free gingiva by precipitating tissular and serum proteins.

INDICATION

- gingival exudation and gingival haemorrhages,
- gingival retraction by in situ impregnation of a cord.

INSTRUCTIONS FOR USE

- Gingival exudation and gingival haemorrhage

Instil **one or two drops** of solution onto the gingiva using a pipette or a cotton swab. Wait 1 to 2 minutes and thoroughly rinse the application site.

- Gingival retraction

First insert a (non-medicated) gingival retraction cord into the sulcus. Carefully impregnate the cord in situ with a few drops of solution, using a pipette. Leave the cord in place for 3 to 5 minutes. Remove it, making sure that no fragment remains into the gingival sulcus. Wash thoroughly.

CONTRAINDICATIONS

Possible allergy to aluminium.

Any history of hypersensitivity or photosensitivity to quinoleic ingredients.

WARNINGS AND PRECAUTIONS FOR USE

- Instil only the number of drops required to obtain the desired result.
- Reduce to a minimum the contact time of the solution with the mucosa.
- As soon as local haemostasis is obtained or the impression performed, thoroughly wash the application site so that there is no potentially irritating aluminium salt left.
- Not to be used for over-impregnation of a medicated cord.
- Do not swallow.

Flammable.

Xi - Irritating to eyes and skin.

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

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