

Product Group Hypodermic Syringes for single use:

Omnifix® Luer Solo, Omnifix® Luer Duo,
Omnifix® Luer Lock Solo, Omnifix® Luer Lock Duo,
Levograf Luer Lock Solo, Levovist Luer Lock,

Product Group Fine-Dosage Syringes for single use:

Omnifix®-F Luer Solo, Omnifix®-F Luer Duo,
Omnifix®-H Luer Solo, AS Luer Solo 1 ml

Customized kits**2.2.1 Description of the device materials (components, raw material)****Syringes**

Device component		Polymer ¹	Polymeric material
1	Barrel	PP	Polypropylene
2	Plunger		
2.1	Plunger 1 ml	PS	Polystyrene
2.2	Plunger 2-50 ml	PP	Polypropylene
3	Piston seal	IR	Polyisoprene
4	External lubricant	--	Silicone fluid

Detached hypodermic needles

see Technical Documentation 39.05.220

¹ Polymer abbreviations according to DIN EN ISO 1043-1 and ISO 1629